
# elements

## *body*

**SN**

.....

BRUKSANVISNING - BRUGERINSTRUKTION - USER INSTRUCTION

BRUKERVEJLEDNING - KÄYTTÖOHJEET

allard

# SVENSK

## Allmän information

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# ENGLISH

## General information

Your Elements Body Orthosis has been prescribed by your clinician to help improve your daily living. It should fit snugly against the body and provide pressure on the trunk and extremities. It is intended to improve your posture and increase stability. For example the hip and shoulder position can be improved. Wearing this Orthosis may also make your speech clearer. The changes may not be instant, you may only see improvement over a period of time.

The fabrics used in Elements Body are carefully selected and are produced & certified according to Oeko-Tex standard 100, product class II. Due to early use of Elements Body in some cases, additional testing has been carried out to meet safety requirements for children under three (3) years. These tests are quite extensive and requirements are high so therefore, only beige fabrics have been tested and approved for infants.

## Purpose of Elements Body

This Orthosis is custom made for you and should not be used by anyone else. It is intended to:

- Increase proprioceptive feedback (awareness of body in space)
- Improve postural control
- Improve body awareness
- Improve fine motor skills
- Improve body posture
- Extremities movement can also be improved.

Other \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

For optimum results it is important that the Orthosis is used every day during activities.

Please follow the instructions provided by your clinician. During rest, the Orthosis will have no function.

Once you have been using the Orthosis for a couple of months you will also notice that the effect of wearing the Orthosis will remain with you even when you are not wearing the suit.

## Advise for applying and wearing Elements Body:

### Applying the Full Suit:

- Open all zippers and closures.
- Pull each leg above the knees one at a time.
- Check the knee panels are centred.
- Pull the suit over the buttocks and hips.
- Insert arms one at a time and then pull the suit over the shoulders.
- Check the elbow panels are centred.
- Check all seams are straight on the arms, legs and sides of the body.
- Pull up each zip one at a time, in any order that you feel most comfortable.

### Wearing the Full Suit:

- By applying a plastic bag over the arm and leg it can simplify donning the Orthosis.
- If you notice any skin irritation or experience a raised temperature, remove the Orthosis and contact your clinician.
- Be observant for signs of discomfort, or change in body temperature. This could be an indication of circulation problems. Immediately remove the Orthosis and contact your clinician.
- The Orthosis should fit snugly. It is normal to see red marks around the neck and arm pits after removing the Orthosis. This should last for about 30-40 minutes. If the red marks do not go away, then contact your clinician
- If the orthosis is overused e.g. more than 10 hours consecutively, you may feel a certain degree of over stimulation that can lead to excessive tiredness.

## Special Instructions for the use of Elements Body for children:

- Be observant that the child does not lose their balance and fall during donning.
- Check the Orthosis is correctly applied and that all reinforcements and zippers are in the correct position.
- When the child grows and / or gains weight the Orthosis may cause more red marks or chaffing on the skin than normal. Contact your clinician to check the fit of the Orthosis.
- Check the size of the Orthosis every 6 months.

# elements body

## Wearing Schedule

- On the first day only use the Orthosis for a maximum of 2 hours at one time.
- If after 2 hours there are no adverse reactions, gradually increase the use by 1 hour every day until you reach the maximum recommended period of use. It is advised that the suit is worn a minimum of 6 hours and a maximum of 10 hours per day. If the suit is worn more than 10 consecutive hours, the patient may in certain conditions feel an overstimulation which may cause tiredness.
- The Orthosis should not be used during the night unless the clinician has instructed this.

## Washing Instructions

Please use a non-biological mild detergent (a synthetic washing powder that does not contain enzymes).

We recommend that Elements Body is washed gently using a laundry bag. Excessive washing may deteriorate the fabric and function of the garment.

Machine wash 40°  
gentle cycle

Dry flat

Do not iron

Do not tumble dry

Do not bleach

Do not Dry Clean

# Support for Better Life!

Camp Scandinavia AB, Karbingatan 38  
SE-254 67 Helsingborg, Sweden  
Phone: +46 42 25 27 01

allard|  
INT'L

[www.allardint.com](http://www.allardint.com)

Tel +46 42 25 27 00  
Fax +46 42 25 27 25

**ALLARD INTERNATIONAL**  
Camp Scandinavia AB  
Karbingatan 38  
SE-254 67 Helsingborg SWEDEN  
[info@allardint.com](mailto:info@allardint.com)