
.....

BRUKSANVISNING - BRUGERINSTRUKTION - USER INSTRUCTION

BRUKERVEJLEDNING - KÄYTTÖOHJEET

allard|_N_T

SVENSK

Allmän information

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ENGLISH

General information

Your HyMo Orthosis has been prescribed by your clinician to help improve your daily living. It should fit snugly against the body and provide pressure on the trunk and extremities. It is intended to improve your posture and increase your joint stability. For example, to improve your hip and shoulder position.

The changes may not be instant, you may only feel slight changes initially. This feeling is a foundation to build on and hopefully you will notice an improvement over the coming months

The fabrics used in HyMo orthosis are carefully selected and are produced & certified according to Oeko-Tex standard 100, product class II.

Intended use

The HyMo orthosis is a textile orthosis in the form of a body and / or trousers with short or long sleeves and legs. The orthosis is intended for patients with Ehlers Danlos Syndrome. The garments are custom-made with a focus on stabilising and limiting the range of motion in the joints.

Indications

Ehlers Danlos Syndrom.

Contraindications

Edema, dislocated joints (not reduced).

Limitations

When used for pool training, a separate set of orthoses is recommended as it is important that orthoses, when used outside the pool, are clean and dry. In order to achieve the best benefit with your orthosis, it is important that it is used every day at the beginning of treatment. Follow the schedule for use in this information or as recommended by your clinician.

Advise for applying and wearing HYMO:

Body

- Open all zippers and the closure at the crotch.
- Step into the garment and pull up over the bottom and the hips.
- Thread in one arm at a time and then pull up over the shoulders. Make sure the elbow panels are centered.
- Pull up the front (long) zipper first and then the rear (short) to get the shoulders in a good position.

Trousers

- Open all zippers.
- Step into the garment and pull each leg separately over your knees.
- Check that the knee panels are centered.
- Pull up the trousers to the waist and close with zipper and Velcro.
- Make sure the seams run straight along the sides of the legs.

Additional advice:

- If you develop a rash or fever, remove the orthosis and contact your clinician.
- Be aware of your skin colour or body temperature change. This may be a sign of circulation problems. Immediately remove the orthosis and contact your clinician.
- The Orthosis should fit snugly. It is normal to see red marks around the neck and arm pits after removing the Orthosis. This should last for about 30-40 minutes. If the red marks do not go away, then contact your clinician.

Warnings:

- Should always be used dry, besides use for pool training.
- Only for single patients use
- Not flame retardant
- Check garments and skin daily before use.

Wearing Schedule

- On the first day only use the Orthosis for a maximum of 2 hours at one time.
- If after 2 hours there are no adverse reactions, gradually increase the use by 1 hour every day.
- Continue to increase usage by one hour per day until you have reached a usage of 8-10 hours a day.
- If you want to use it longer, take a break for 2 hours. If you use it during the night, take a break for 2 hours.

Washing Instructions

- Please use a perfume free mild detergent.
- We recommend that the HyMo orthosis is washed gently using a laundry bag. Close all zippers and Velcro straps before washing.
- Wash with similar colours.
- Do not use fabric conditioner
- Do not dry in direct sunlight
- Do not dry over the radiator

Machine wash 40°

Do not iron

Do not tumble dry

Do not bleach

Support for Better Life!

Camp Scandinavia AB, Karbingatan 38
SE-254 67 Helsingborg, Sweden
Phone: +46 42 25 27 01

©Camp Scandinavia April 2019

SVERIGE

042-25 27 01
042-25 27 25
info@camp.se

DANMARK

43 96 66 99
43 43 22 66
info@camp.dk

FINLAND

09-350 76 30
09-350 76 338
info@camp.fi

NORGE

23 23 31 20
23 23 31 21
info@camp.no

UK & IRELAND

+353-42 932 8177
+353 42 932 8182
info@allarduk.co.uk

USA

888-678-6548
800-289-0809
info@allardusa.com

INTERNATIONAL

+46 42 25 27 01
+46 42 25 27 25
info@allardint.com

Support for Better Life! ◦ www.allardsupport.com