ENGLISH

INDICATIONS:

- Mild to moderate unicompartmental knee osteoarthritis with ligament instabilities
- Increased rehabilitative and functional support for diagnoses that require reduced medial or lateral joint loads

WARNINGS:

If you experience pain, swelling, sensation changes, or unusual reactions while using this product, contact VQ OrthoCare's Patient Care department at 800.452.7993 or contact the provider that delivered your brace or consult a physician.

Any serious incident that has occurred in relation to the device should be reported to VQ OrthoCare and the competent authority of the Member State in which the user and/or patient is established.

FOR SINGLE PATIENT USE ONLY.

NOTICE: While every effort has been made using

state-of-the-art techniques to obtain maximum compatibility of function, strength, durability and comfort, this device is only one element in the overall treatment program administered by a medical professional. There is no guarantee that injury will be prevented through use of this device.

CAUTION: Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of this device.

WARRANTY: VQ OrthoCare warrants all its products from the original date of purchase against defects in materials and workmanship. Normal wear and tear during use of a product is not considered a defect. Contact VQ OrthoCare at 949.261.3000 for specific product warranty information. Outside the US, please contact the provider that delivered your brace.

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EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands


Subject to one or more patents listed on
www.vqorthocare.com/patents

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








4a

ENGLISH

Please read these instructions in their entirety before applying the brace.

STEP 1. POSITION BRACE

- Open the brace by undoing the pushbuttons and separating the front Velcro closures.
- While sitting at edge of chair, bend knee slightly and apply brace to lateral side (outside) of leg.
- NOTE:** The brace must be applied directly to the skin.
- Align center of hinge with top of kneecap (1a) and position hinge slightly behind midline on the outside of the leg (1b).

NOTE FOR INITIAL FITTING: After above positioning, the brace should be adjusted to neutral. This is done by setting the brace angle (see step 4 below) so the frame and hinge rest on leg without gapping or pressure (1c). Medial or lateral compartment correction will be adjusted in step 4.

STEP 2. FASTEN WRAP

- Secure the top and then the bottom wrap Velcro closures (2a) being sure to maintain hinge position. Use finger pockets on wrap to aid grasp.

STEP 3. FASTEN STRAPS

Medial Compartment OA

- Pull the **TOP** strap behind the thigh and insert pushbutton into **TOP** buckle (3a). Tighten strap and secure Velcro pull-tab (3b). Do not overtighten.
- Repeat above step for **BOTTOM** strap securing it to the **BOTTOM** buckle (3c). This strap should be tightened very securely, but comfortably, as it keeps the brace in position.
- View from behind the knee (3d).

Lateral Compartment OA (Criss-Cross Strapping)

- Pull the **TOP** strap downward behind the knee and insert pushbutton into **BOTTOM** buckle. Attach the Velcro pull-tab to the strap keeping strap slightly loose.
- Repeat above step for **BOTTOM** strap securing it to the **TOP** buckle.
- View from behind the knee (3e).

Step 4. Apply Correction

- With the brace in neutral (from step 1 initial fitting), bend the knee to 90° (4a).

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- For medial compartment relief, use the adjustment tool to turn the Q-hinge screw (4a) 1/4 turn counterclockwise to move the hinge in toward the knee.
- For lateral compartment relief, use the adjustment tool to turn the Q-hinge screw (4a) 1/4 turn clockwise to move the hinge away from the knee. Important note: Adjust the straps to remove slack but remain slightly loose. Over-tightening may cause excessive pressure upon knee extension.
- The above is a good starting point for initial use.
- For initial fitting and adjustments, stand and walk a short distance to ensure comfort and proper correction after each adjustment. If there is any discomfort in the knee joint or leg, turn the angle adjustment back slightly.
- If more relief is required, turn Q-hinge screw no more than an additional 1/4 turn. Walk after the adjustment and check as above.
- Strap readjustment may be needed after Q-hinge adjustments.
- Over-adjusting the Q-hinge may cause discomfort, rotation or skin pressure issues. Adjust conservatively.

NOTE: After initial application and adjustment, the straps do not need to be loosened and tightened for each application. The brace can be applied and removed using just the buckles.

“WEAR TIME”

Gradually increasing brace wear time will allow your body to adjust to the brace. Example: Start with 30 - 60 minutes on day one and double the amount of time each day. Pain caused by osteoarthritis may take time to subside. It is important that you routinely wear the brace when you're on your feet as it will reduce joint loading and allow you to begin feeling relief.

CARE AND CLEANING

Gently hand wash the brace in cold water with mild detergent. Rinse well and air dry. Do **NOT** machine wash or dry.

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© 2022 VisionQuest Industries Inc.