BOSTON RYGGKORSETTER

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BOSTON SPINAL ORTHOSIS

- Wear & Care Guide (ENG)

Your doctor/orthotist has prescribed this brace to assist with a particular condition. As the treatment can vary from case to case your doctor or orthotist will provide specific information about:

- Applying your Orthosis
- Weaning (Acclimatization) schedule to achieve optimum hours of daily wearing.
- Duration of use (optimum wearing schedule)
- Acceptable activities
- Activities to avoid
- Who to contact in case of problem

Cleaning your Orthosis:

- It is important to clean the orthosis on a daily basis. For best results use rubbing alcohol. Alternatively use a soapy cloth and rinse with a damp cloth and towel dry.
- Always check your orthosis before putting it back on to ensure that it undamaged and free from dirt.

Skin Care

It is **very important** to prevent skin breakdown (sore, red, raw skin). The skin underneath the orthosis needs to be toughened up, especially where the orthosis applies the most pressure.

- Bathe on a daily basis.
- Pay special attention to pink areas of the skin where orthosis pressure is highest.
- Always wear a snug-fitting t-shirt (without seams) under your orthosis and ensure to stretch the fabric to ensure that it is free from wrinkles and also that it is dry before putting on your brace (available from your orthotist).
- If skin irritation/breakdown persists, contact your fitter or clinic for advise.

NOTE:

Any special Doctor's / Orthotist orders should always supersede the general guideline information given in this guide.

This product is to be fitted by an appropriately trained and qualified professional and with reference to the specific product Technical Manual, available to view and download at www.allarduk.co.uk