allard

Support for better life!

User Guide

Custom Scoliosis Night Braces

Intended Use

The Boston Scoliosis Night Brace and the Hybrid Night Brace are custom spinal orthoses used at night-time in the non-operative and conservative treatment of Idiopathic Scoliosis which primarily used for adolescents.

Indication / Contraindication

Boston Scoliosis Night Brace	Hybrid Night Brace
Indications:	Indications:
Adolescent Idiopathic Scoliosis (AIS)	Curves with Cobb angle between 20° and 45° both single and double. Risser 0-4 Procrastination for surgery Apex Thoracic no higher than TH 6-7 Curves below 25° where the gibbus is greater than 7° and familial Scoliosis is known.
Contraindications:	Contraindications:
Thoracic Lordosis Severe obesity	None known

General Safety Instruction

- This orthosis is for multiple use by a single patient
- Only appropriately trained and qualified Clinicians should fit your orthosis with specific reference to the Technical Manual available to view and download at www.allarduk.co.uk
- Patients should contact their Clinician in event of any persistent skin irritation
- In the event of incorrect use product liability is excluded

Application Information

Your Clinicians instruction will always supersede the information given in this instruction!

How to put on your Boston Scoliosis Night Brace:

The orthosis is more effective when put on in a laying down position as it is easier to tighten the straps!

Put on the orthosis while standing up and hold it against the body. Then lay on your back. Ensure that the waist (the pads located at the centre on the inside of the orthosis) fits above the hip bone and under the ribs.

Begin to tighten the middle strap on the front of the orthosis. If there are only two straps begin with the lower strap. Then tighten the top strap and finally the lower. Adjust all straps so they are evenly tightened. Always pull the straps as tight as possible to avoid the orthosis sliding and reduce the risk of skin irritation.

How to put on your Hybrid Night Brace:

Put on the orthosis while standing up and hold it against the body. Then lay on your back. Ensure that the waist (the pads located at the centre on the inside of the orthosis) fits above the hip bone and under the ribs.

Begin to tighten the top strap on the front of the orthosis. Then tighten the bottom strap and final the middle strap. If there are only two straps begin with the top strap and then tighten the bottom strap. Adjust all straps so they are evenly tightened. Always pull the straps as tight as possible to avoid the orthosis sliding and reduce the risk of skin irritation.

Getting used to your Custom Night Brace:

It can take up to a week to get used to wearing the orthosis. The first few nights you may wake up several times during the night. If this happens remove the orthosis and sleep the rest of the night without the orthosis on. Repeat this process each night and try to build up wearing time so that you eventually sleep through the whole night wearing the orthosis. If several weeks past and you are still not able to sleep through the night wearing the orthosis, please contact the clinician who fitted the orthosis.

Side Effects

The skin under the orthosis should be checked regularly for irritation. Due to the abnormal pressure of the orthosis on your skin while wearing for extended lengths of time your skin will need time to adjust.

We recommend that you:

- Shower or bathe daily;
- Pay attention to where the skin has become pink (where the most pressure is);
- Wear a seamless, tight-fitting vest under the orthosis which covers the body under the orthosis.

Contact your clinician if any irritation that occurs does not heal.

Maintenance and Cleaning

It is important to clean the orthosis daily. Use a damp cloth and soapy water to wipe it clean. Dry thoroughly. For best results, clean the orthosis in the morning and allow to dry during the day.

CE Conformity

The product satisfies the requirement regulation (EU) 2017/745 of the European Parliament and of the Council (MDR) and does not bear a CE mark being a custom made orthosis. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

Warranty and Guarantee

We recommend that the orthosis is renewed / changed if it no longer fits or is damaged.

Warranty is provided under the terms of sale and supply of Allard Support UK Ltd provided that the above conditions are met.

Storage and Disposal

When not in use store your orthosis safely to avoid accident damage. Refer to your local disposal regulations with reference to the material content as per the material content symbol on your orthosis

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