

OACTIVE® 2

OActive® 2 Osteoarthritis Knee Brace is a lightweight, and low-profile upright osteoarthritis knee brace with unique varus / valgus angle adjustment to optimize pain relief. Laterally applied frame and hinge corrects for either medial or lateral compartment and avoids contralateral limb contact during gait.

The unique Q-joint is adjusted clockwise or counterclockwise to achieve the desired relief. Start with the orthosis in the neutral position, with the knee bent at 90°, and turn 1/8 turn at a time to adjust the relief. OActive® 2 can also be set to limit the range of motion for extension and flexion independently. The orthosis can also be set for free movement. Comes with 5° stop installed.

Available in Set-Fit or Push Button Buckle System in either:

- Standard hinge for patients weighing under 136Kgs
- Heavy Duty hinge for patients weighing over 136Kgs

Features & Benefits

- Lightweight, low-profile, single upright OA brace
- Unique Set-Fit Buckle System eliminates the hassle of Velcro® adjustments each time the brace is applied and removed.
- Single hand Push Button Buckle for patients with limited hand function
- Maximized 3 points of leverage for medial compartment
- Laterally applied frame corrects for either medial or lateral compartment
- Lateral application avoids opposite limb during gait
- Supplied with Suspension Wrap or Tibia Cuff

Recommended Range Of Application

Mild to severe unicompartmental osteoarthritis.



Item No.	Model	Size	L/R
550001000	Standard Hinge Push Button with Suspension Wrap	<136Kg	Left
550002000	Standard Hinge Push Button with Suspension Wrap	<136Kg	Right
550171000	Heavy Duty Hinge Push Button with Suspension Wrap	>136Kg	Left
550172000	Heavy Duty Hinge Push Button with Suspension Wrap	>136Kg	Right
550291000	Standard Hinge Push button with Tibia Cuff	<136Kg	Left
550292000	Standard Hinge Push button with Tibia Cuff	<136Kg	Right