

BlueROCKER® 2½

The Foot Drop Orthosis BlueROCKER® 2 ½ has the same dynamic functions and stability as a BlueROCKER® 2.0 with the exception that the 2 ½ has a lower heel height, 7 mm, to adapt to modern shoes. The toe lift is also slightly lower, which results in more space in the toe box. Since the start in 1997 the Allard AFO has been the market leader with its unique functions and dynamic properties. BlueROCKER® 2 ½ has all MikroFIX pre-applied to simply fitment of the padding. The strap solution is the same as on 2.0, made from a durable and strong material. The straps can be applied for either inside or outside closure on the orthosis. The foot plate on BlueROCKER® 2 ½ gives more space in the toe box. It is always important to adapt any Allard AFO to the heel height of the shoe. Read more about this in the Allard AFO Professional Instructions.

BLUEROCKER® 2 ½



Recommended Range Of Application

BlueROCKER® is designed to manage footdrop in conditions such as Stroke, Multiple Sclerosis, Post Polio Syndrome, Muscular Dystrophy, Spinal Cord Injuries, Traumatic Brain Injuries, Guillian-Barre Syndrome, Charcot-Marie-Tooth, Myelomeningocele, Neuropathy or Cerebral Palsy. BlueROCKER® 2.0 is designed to support gait in conditions such as Posterior Tibialis Tendon Dysfunction (PTTD) and toe amputations. BlueROCKER® can also be used for partial foot amputations, most proximal level is Chopart.

Contraindications

BlueROCKER® 2 ½ should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits (e.g. quadriceps spasticity, genu valgum or varum, genu recurvatum), severe spasticity.

Item No.	Model	Size	L/R	Footplate Length	Height
289701010	BlueROCKER® 2½	X-Small	Left	21 cm	36 cm
289701011	BlueROCKER® 2½	Small	Left	23 cm	38 cm
289701012	BlueROCKER® 2½	Medium	Left	24,5 cm	40,5 cm
289701013	BlueROCKER® 2½	Large	Left	27 cm	43 cm
289701014	BlueROCKER® 2½	X-Large	Left	28,5 cm	45 cm
289702010	BlueROCKER® 2½	X-Small	Right	21 cm	36 cm
289702011	BlueROCKER® 2½	Small	Right	23 cm	38 cm
289702012	BlueROCKER® 2½	Medium	Right	24,5 cm	40,5 cm
289702013	BlueROCKER® 2½	Large	Right	27 cm	43 cm
289702014	BlueROCKER® 2½	X-Large	Right	28,5 cm	45 cm