

Suit Reinforcements 39310, 39311, 39312

elements
body



SP
Standard Panel

To assist hip abduction and back extension

LP-X
Lumbar Panel Extended

To assist postural control

HPL
Hip Panel (left side)

To assist hip control and abduction

HPR
Hip Panel (right side)

To assist hip control and abduction

LP
Lumbar Panel

To assist postural control

MLL
Medial Leg left

To assist adduction and internal rotation

LTPL
Lateral Trunk Panel

To assist correction of LEFT lateral flexion

MLR
Medial Leg Right

To assist adduction and internal rotation


LTPR
Lateral Trunk Panel

To assist correction of RIGHT lateral flexion

Three reinforcements are included in the price. If the box for Standard Panel is ticked this will count as one, and another two can be added.

Suit Reinforcements 39310, 39311, 39312

elements
body



PHE-L

Posterior Hip Extension Left


To assist hip extension on LEFT side

PHE-R

Posterior Hip Extension Right

To assist hip extension on RIGHT side

Long Legs and $\frac{3}{4}$ Legs only



AKE-L

Anterior Knee Extension Left


To assist knee extension on LEFT side

AKE-R

Anterior Knee Extension Right

To assist knee extension on RIGHT side

Long Legs and $\frac{3}{4}$ Legs only



PKF-L

Posterior Knee Flexion Left

To assist knee flexion on LEFT side

PKF-R






Posterior Knee Flexion Right

To assist knee flexion on RIGHT side

Long Legs and $\frac{3}{4}$ Legs only



Suit Reinforcements 39310, 39311, 39312

elements
body



Suit Reinforcements 39310, 39311, 39312

elements
body




IRPL-R

Long Internal Rotation Panel

To assist internal rotation at
RIGHT hip

Only for Long Legs




IRPL-L

Long Internal Rotation Panel

To assist internal rotation at
LEFT hip

Only for Long Legs




ERPL-R

Long External Rotation Panel

To assist external rotation at
RIGHT hip

Only for Long Legs



ERPL-L

Long External Rotation Panel


To assist external rotation at
LEFT hip

Only for Long Legs

Suit Reinforcements

39310, 39311, 39312

elements
body



PPP
Posterior Pelvic Panels
To assist anterior tilt

APP
Anterior Pelvic Panels
To assist posterior tilt

SBEL
Shoulder Back Extension Left
To assist retraction of the scapulae on the left side and extension of the back on the RIGHT side

SBER
Shoulder Back Extension Right
To assist retraction of the scapulae on the right side and extension of the back on the LEFT side

SPL
Scoliosis Panels Left
To assist correction of right thoracic and LEFT lumbar curves

CENTER
Right thoracic
Left Lumbar

Please ensure that the correct reinforcement is selected.

PBP
Posterior Back Panel
To assist extension of the trunk

Model: Paediatric, Male
Model: Female

SPR
Scoliosis Panels Right
To assist correction of left thoracic and RIGHT lumbar curves

CENTER
Left thoracic
Right Lumbar


Please ensure that the correct reinforcement is selected.

AAP
Anterior Abdominal Panels
To assist posterior tilt

Three reinforcements are included in the price. If the box for Standard Panel is ticked this will count as one, and another two can be added.

39310, 39311, 39312 Suit Reinforcements

elements
body



SPL-C3
Scoliosis Panels Left

To assist correction of
LEFT C-curve and right
trunk rotation

CENTER

SPR-C3
Scoliosis Panels Right

To assist correction of
RIGHT C-curve and left
trunk rotation

CENTER

SPL-C2
Scoliosis Panels Left

To assist correction
of LEFT C-Curve and
manage RIGHT Lateral
trunk flexion

CENTER

SPR-C2
Scoliosis Panels Right

To assist correction of
RIGHT C-Curve and ma-
nage LEFT Lateral trunk
flexion

CENTER

SPL-C1
Scoliosis Panels Left


To assist correction of
LEFT C-curve with mild
right trunk rotation

CENTER

SPR-C1
Scoliosis Panels Right

To assist correction of
RIGHT C-curve with
mild left trunk rotation

CENTER



AS-L
Anterior
Supination Left

To resist supination
on left side

Long Sleeve Only



AS-R
Anterior
Supination Right

To resist supination
on right side


Long Sleeve Only



PP-L
Posterior
Supination Left

To resist pronation
on left side

Long Sleeve Only



PP-R
Posterior
Pronation Right


To resist pronation
on right side

Long Sleeve Only

Three reinforcements are included in the price. If the box for Standard Panel is ticked this will count as one. Another two can be added. If Standard Panel is not ticked, three reinforcements can be chosen.

Suit Reinforcements 39312


elements
body




AEP left
Anterior Elbow
Panel Left
To resist elbow hyper-
extension on left side
Long Sleeve Only



AEP right
Anterior Elbow
Panel Right
To resist elbow hyper-
extension on right side
Long Sleeve Only



PEP left
Posterior Elbow
Panel Left
To resist elbow flexion
on left side
Long Sleeve Only



PEP right
Posterior Elbow
Panel Right
To resist elbow flexion
on right side
Long Sleeve Only

Three reinforcements are included in the price. If the box for Standard Panel is ticked this will count as one. Another two can be added. If Standard Panel is not ticked, three reinforcements can be chosen.

allard

Allard Support UK Ltd.

Coes Road, Dundalk, Co. Louth, A91 E868 | +353 42 932 8177 | info@allarduk.co.uk | allarduk.co.uk