



elements *body*

TRIAL SLEEVE

The Trial sleeve has been designed to be used in the assessment process for Elements Body MTM upper limb garments. It enables the clinician to evaluate the users response to the application of a dynamic textile garment in preparation for a made to measure upper limb garment. The Trial sleeve is only to be used for an assessment period of 2-4 weeks to assist the clinician in identifying the clinical effectiveness and tolerance prior to the provision of an Elements Body MTM Garment.



Item No.	Size	L/R	Wrist cm	Elbow cm	Biceps cm
39344 0008	XXX-Small	Bilat.	11-12	16-18	16-18
39344 0009	XX-Small	Bilat.	12-13	18-20	18-20
39344 0010	X-Small	Bilat.	13-14	20-22	20-22,5
39344 0011	Small	Bilat.	14-15	21-23	22,5-25
39344 0012	Medium	Bilat.	14-16	22-25	24-27
39344 0013	Large	Bilat.	16-18	25-28	27-30
39344 0014	X-Large	Bilat.	18-20	28-31	30-33
39344 0015	XX-Large	Bilat.	20-22	31-34	33-37