


39330, 39331, 39335, 39336 Glove
39340 Sleeve
Reinforcements

elements
body



AS-L
Anterior Supination Left

*To resist supination
on left side*

The diagram shows two human figures from the waist up. The left figure has orange and white diagonal stripes on the front of its left forearm. The right figure has the same stripes on the front of its right forearm.

AS-R
Anterior Supination Right

*To resist supination
on right side*

The diagram shows two human figures from the waist up. The left figure has orange and white diagonal stripes on the front of its right forearm. The right figure has the same stripes on the front of its left forearm.

PP-L
Posterior Pronation Left

*To resist pronation
on left side*

The diagram shows two human figures from the waist up. The left figure has orange and white diagonal stripes on the back of its left forearm. The right figure has the same stripes on the back of its right forearm.

PP-R
Posterior Pronation Right

*To resist pronation
on right side*



The diagram shows two human figures from the waist up. The left figure has orange and white diagonal stripes on the back of its right forearm. The right figure has the same stripes on the back of its left forearm.

Three reinforcements are included in the price. If the box for DW is ticked this will count as one. Another two can be added. If DW is not ticked, three reinforcements can be chosen.

allardku

39330, 39331, 39335, 39336 Glove
39340 Sleeve
Reinforcements


elements
body



AEP-L

Anterior Elbow Panel Left


To resist elbow hyperextension on left side



AEP-R

Anterior Elbow Panel Right


To resist elbow hyperextension on right side



PEP-L

Posterior Elbow Panel Left

To resist elbow flexion on left side



PEP-R

Posterior Elbow Panel Right



To resist elbow flexion on right side

Three reinforcements are included in the price. If the box for DW is ticked this will count as one. Another two can be added. If DW is not ticked, three reinforcements can be chosen.

allard|ku

39330, 39331, 39335, 39336 Glove
39340 Sleeve
Reinforcements


elements
body




VW
Volar Wrist
To resist wrist extension




VWx
Volar Wrist
 extended to end of
 Fingers
To resist wrist extension




DW
Dorsal Wrist
To resist Wrist flexion



DWx
Dorsal Wrist
 extended to end of
 Fingers
To resist Wrist flexion



TP
Thumb Panel
To assist extension and abduction of the Thumb



DP
Deviation Panel
To correct Ulnar Deviation

Three reinforcements are included in the price.
 If the box for DW is ticked this will count as one, and another two can be added.
 If DW is not ticked, three other reinforcements can be chosen.