


**39300, 39301 Vest
39305, 39306 Body
Reinforcements**

elements
body



SP
Standard Panel

To assist hip external rotation and back extension

LTPR
Lateral Trunk Panel

To assist correction of RIGHT lateral flexion

LP-X
Lumbar Panel Extended

To assist postural control

SPL
Scoliosis Panels Left

To assist correction of right thoracic and LEFT lumbar curves

Please ensure that the correct reinforcement is selected.

LP
Lumbar Panel

To assist postural control

SPR
Scoliosis Panels Right

To assist correction of left thoracic and RIGHT lumbar curves


Please ensure that the correct reinforcement is selected.

LTPL
Lateral Trunk Panel

To assist correction of LEFT lateral flexion

39300, 39301 Vest
39305, 39306 Body
Reinforcements

elements
body



SPL-C3
 Scoliosis Panels Left
 To assist correction of *LEFT* C-curve and right trunk rotation

SPR-C3
 Scoliosis Panels Right
 To assist correction of *RIGHT* C-curve and left trunk rotation

SPL-C2
 Scoliosis Panels Left
 To assist correction of *LEFT* C-Curve and manage *RIGHT* Lateral trunk flexion



SPR-C2
 Scoliosis Panels Right
 To assist correction of *RIGHT* C-Curve and manage *LEFT* Lateral trunk flexion

SPL-C1
 Scoliosis Panels Left
 To assist correction of *LEFT* C-curve with mild right trunk rotation

SPR-C1
 Scoliosis Panels Right
 To assist correction of *RIGHT* C-curve with mild left trunk rotation

39300, 39301 Vest 39305, 39306 Body Reinforcements

elements
body



SBEL


Shoulder Back Extension Left
To assist retraction of the left scapulae and back extension on the RIGHT side

SBER

Shoulder Back Extension Right
To assist retraction of the right scapulae and back extension on the LEFT side

SBEL

SBER




PBP

Posterior Back Panel

To assist trunk extension

Model:
Paediatric,
Male

Model:
Female



AAP

Anterior Abdominal Panel



To assist trunk flexion

For Body Only

Three reinforcements are included in the price. If the box for Standard Panel is ticked this will count as one. Another two can be added. If Standard Panel is not ticked, three reinforcements can be chosen.

3930I Vest
39306 Body
Reinforcements


elements
body



AEP-L
Anterior Elbow Panel Left

To resist elbow hyperextension on left side


Long Sleeves Only



AEP-R
Anterior Elbow Panel Right

To resist elbow hyperextension on right side


Long Sleeves Only



PEP-L
Posterior Elbow Panel Left

To resist elbow flexion on left side

Long Sleeves Only



PEP-R
Posterior Elbow Panel Right


To resist elbow flexion on right side

Long Sleeves Only

Three reinforcements are included in the price. If the box for Standard Panel is ticked this will count as one. Another two can be added. If Standard Panel is not ticked, three reinforcements can be chosen.

3930I Vest 39306 Body Reinforcements

elements
body



AS-L
Anterior Supination Left

*To resist supination
on left side*


Long Sleeves Only



AS-R
Anterior Supination Right

*To resist supination
on right side*


Long Sleeves Only



PP-L
Posterior Pronation Left

*To resist pronation
on left side*

Long Sleeves Only



PP-R
Posterior Pronation Right

*To resist pronation
on right side*

Long Sleeves Only

Three reinforcements are included in the price. If the box for Standard Panel is ticked this will count as one. Another two can be added. If Standard Panel is not ticked, three reinforcements can be chosen.

allard|KC