

HyMo – Dynamic Textile Orthosis for Hypermobility

HyMo is a made-to-measure textile orthosis developed to provide external support for hypermobile joints and surrounding soft tissues, particularly in people with Ehlers-Danlos syndrome (EDS).

It fits closely to the body and can be customised with targeted reinforcement to help improve alignment, enhance joint stability, and reduce pain, while supporting function in everyday activities.

Product options

HyMo garments are available as a Body and/or Trousers, with short, three-quarter or long sleeves and legs.

Recommended Range Of Application

Hypermobile Ehlers-Danlos syndrome (hEDS) and Hypermobility Spectrum Disorders (HSD), characterised by generalised joint hypermobility and associated musculoskeletal symptoms.



Material & function

HyMo is intended for individuals with joint hypermobility affecting joints and associated soft tissues, including those diagnosed with Ehlers-Danlos syndrome (EDS) or Hypermobility Spectrum Disorders (HSD). It is provided exclusively as a made-to-measure orthosis and is designed to support joint stability and help control excessive range of motion.

The HyMo material contains a high proportion of Lycra®, providing stretch and recovery while maintaining dynamic support. This enables functional movement during strengthening and conditioning activities, including hydrotherapy. The fabric is breathable and quick-drying, supporting comfort during daily use.

HyMo is classified as a Class I Medical Device and is produced and certified in accordance with MDD 93/42/EEC and MDR requirements.

The material is certified according to Oeko-Tex® Standard 100, Class II.

Prescription

HyMo garments should only be prescribed by clinicians who have completed the Allard training programme. A range of garment styles and targeted reinforcement options are available to meet individual patient requirements. Order forms, clinical guidelines, and measurement documentation are available to download at allarduk.co.uk.




Further product information is available via the QR code.



“My HyMo-Body helps me function better in my daily life and I feel more stable and have less pain”

Maria, EDS patient

Item No. (Female/Male)	Model	Color
394000000/394010000	Body with short sleeves	Black
394020000/394030000	Body with long sleeves	Black
394080000/394090000	Body with 3/4 sleeves	Black
394040000/394050000	Trousers with short legs	Black
394060000/394070000	Trousers with long legs	Black
394100000/394110000	Trousers with 3/4 legs	Black



allard

Allard Support UK Ltd.

Coes Road, Dundalk, Co. Louth, A91 E868 | +353 42 932 8177 | info@allarduk.co.uk | allarduk.co.uk