

## ASO® EVO™


Stable low-profile ankle support with stabilising front lacing. The rigid nylon straps mimic taping and prevent the foot from supinating or pronating. On the ASO EVO, the nylon straps around the calcaneus provide greater resistance to calcaneal inversion/eversion than the standard ASO.

The upper part is reinforced with a plastic cuff, enhancing protection against high ankle sprains. For added comfort, the orthosis features new reinforced CoolFlex padding along the inside, next to the Achilles tendon.

### Recommended Range Of Application

Acute or chronic ankle instability, preventive use during sports activities, high ankle sprains, and subtalar joint instability.

### Measurement for donning



Item.no	Model	Size	Measurement	L/R	Colour
281890010	ASO® EVO™	X-Small	25-27,5 cm	Bilateral	Black
281890011	ASO® EVO™	Small	27,5-30 cm	Bilateral	Black
281890012	ASO® EVO™	Medium	30-32,5 cm	Bilateral	Black
281890013	ASO® EVO™	Large	32,5-35 cm	Bilateral	Black
281890014	ASO® EVO™	X-Large	35-37,5 cm	Bilateral	Black
281890015	ASO® EVO™	XX-Large	37,5-40 cm	Bilateral	Black