

## Ryno Lacer® II

Highly stable wrist and thumb orthosis with mouldable aluminium splints on the dorsal, volar, and radial sides. The radial thumb splint extends along the full length of the splint and effectively immobilises the thumb. Made from durable suede-like material with a breathable lining.

The lacing locking system is easy to open and close with one hand. Ryno Lacer II is an updated and improved version of the Ryno Lacer, featuring enhanced palm padding, increased thumb cushioning, and an updated Velcro system on the underside of the orthosis.

### Recommended Range Of Application

Post-traumatic, post-operative, rheumatologic conditions, and spasticity.



Item.no	Model	Size	Measurement	L/R	Colour
285761010	Ryno Lacer® II	X-Small	<14 cm	Left	Black
285761011	Ryno Lacer® II	Small	14-16 cm	Left	Black
285761012	Ryno Lacer® II	Medium	16-18 cm	Left	Black
285761013	Ryno Lacer® II	Large	18-20 cm	Left	Black
285761014	Ryno Lacer® II	X-Large	20-22 cm	Left	Black
285762010	Ryno Lacer® II	X-Small	<14 cm	Right	Black
285762011	Ryno Lacer® II	Small	14-16 cm	Right	Black
285762012	Ryno Lacer® II	Medium	16-18 cm	Right	Black
285762013	Ryno Lacer® II	Large	18-20 cm	Right	Black
285762014	Ryno Lacer® II	X-Large	20-22 cm	Right	Black