

## Tripod™ Wrist Lacer II

The Tripod Wrist Lacer II is a supportive wrist brace designed for enhanced immobilization and easy application. Ulnar, palmar, and dorsal aluminum stays minimize wrist flexion-extension and also help control radial-ulnar deviation for effective stabilization. The thumb strap features an additional hook-and-loop interface for improved security, and the ulnar and palmar stays can be removed to support a progressive treatment approach.

A patented lace closure design distributes equal tension across all laces, creating uniform compression for outstanding support. The tension can be easily adjusted for patients with disproportionate forearm and wrist circumferences. Made from a durable polymer felt liner and perforated suede that allows moisture to evaporate, enhancing comfort and patient compliance.

### Recommended Range Of Application

Post-fracture care, post-operative care, carpal tunnel syndrome, wrist fractures, or wrist strain



Item.no	Model	Size	Measurement	L/R	Colour
285811010	Tripod™ Wrist Lacer II	X-Small	12-14 cm	Left	Black
285811011	Tripod™ Wrist Lacer II	Small	14-16 cm	Left	Black
285811012	Tripod™ Wrist Lacer II	Medium	16-18 cm	Left	Black
285811013	Tripod™ Wrist Lacer II	Large	18-20 cm	Left	Black
285811014	Tripod™ Wrist Lacer II	X-Large	20-22 cm	Left	Black
285812010	Tripod™ Wrist Lacer II	X-Small	12-14 cm	Right	Black
285812011	Tripod™ Wrist Lacer II	Small	14-16 cm	Right	Black
285812012	Tripod™ Wrist Lacer II	Medium	16-18 cm	Right	Black
285812013	Tripod™ Wrist Lacer II	Large	18-20 cm	Right	Black
285812014	Tripod™ Wrist Lacer II	X-Large	20-22 cm	Right	Black