

## Wrist Lacer™ II, Short

Highly stable wrist orthosis with mouldable aluminium splints on both the dorsal and volar sides. The splints can be removed and shaped as needed. A cushion on the inside of the orthosis sleeve protects the tendons at the base of the thumb near the distal lacing. Made from durable suede-like material with a breathable lining.

Locking system with a unique design that provides even compression across all four straps. Compression can be easily adjusted for patients with atypical forearm and wrist shapes. One-handed operation of the lacing system. Adjustable thumb grip.

### Recommended Range Of Application

Carpal tunnel syndrome, post-traumatic, and post-operative use.



Item.no	Model	Size	Measurement	L/R	Colour
285781009	Wrist Lacer II, Short	XX-Small	<12 cm	Left	Black
285781010	Wrist Lacer II, Short	X-Small	12-14 cm	Left	Black
285781011	Wrist Lacer II, Short	Small	14-16 cm	Left	Black
285781012	Wrist Lacer II, Short	Medium	16-18 cm	Left	Black
285781013	Wrist Lacer II, Short	Large	18-20 cm	Left	Black
285781014	Wrist Lacer II, Short	X-Large	20-22 cm	Left	Black
285781015	Wrist Lacer II, Short	XX-Large	22-25 cm	Left	Black
285782009	Wrist Lacer II, Short	XX-Small	<12 cm	Right	Black
285782010	Wrist Lacer II, Short	X-Small	12-14 cm	Right	Black
285782011	Wrist Lacer II, Short	Small	14-16 cm	Right	Black
285782012	Wrist Lacer II, Short	Medium	16-18 cm	Right	Black
285782013	Wrist Lacer II, Short	Large	18-20 cm	Right	Black
285782014	Wrist Lacer II, Short	X-Large	20-22 cm	Right	Black
285782015	Wrist Lacer II, Short	XX-Large	22-25 cm	Right	Black