

## Wrist Lacer™ II, Short

Highly stable long wrist orthosis with mouldable aluminium splints on both the dorsal and volar sides. The splints are removable and can be shaped as required. A cushion on the inside of the orthosis sleeve protects the tendons at the base of the thumb near the distal lacing. Made from durable suede-like material with a breathable lining. Unique locking system providing even compression across all four straps. Compression can be easily adjusted for patients with atypical forearm or wrist shapes. One-handed operation of the lacing system. Adjustable thumb support.

### Recommended Range Of Application

Post-traumatic, post-operative, rheumatologic conditions, and spasticity.



Item.no	Model	Size	Measurement	L/R	Colour
285791010	Wrist Lacer II, Long	X-Small	12-14 cm	Left	Black
285791011	Wrist Lacer II, Long	Small	14-16 cm	Left	Black
285791012	Wrist Lacer II, Long	Medium	16-18 cm	Left	Black
285791013	Wrist Lacer II, Long	Large	18-20 cm	Left	Black
285791014	Wrist Lacer II, Long	X-Large	20-22 cm	Left	Black
285792010	Wrist Lacer II, Long	X-Small	12-14 cm	Right	Black
285792011	Wrist Lacer II, Long	Small	14-16 cm	Right	Black
285792012	Wrist Lacer II, Long	Medium	16-18 cm	Right	Black
285792013	Wrist Lacer II, Long	Large	18-20 cm	Right	Black
285792014	Wrist Lacer II, Long	X-Large	20-22 cm	Right	Black