

SELECTION[®]

Textile wrist- and
thumb support



Functional models that cover your needs!

Wrist-orthoses

SELECTION®

SELECTION® Wrist support Soft & Rigid

SELECTION Wrist support is a development from an older very successful product within our product line - Carpus Flexi. These two new models keep the lacing system together with new added material features.

SELECTION wrist supports comes in two versions, one with a thumb support, rigid version, and the other one with a functional low cut in the palm of the hand, soft version. Both have the adjustable lacing on the ulnar side.

Each Selection Wrist support have the new CS-spoon inserted on the palmar side of the orthosis. The CS-spoon is removable.

| Art.nr. | Colour | Model | Left side/ Right side | Size | Wrist circumference |
|-------------|--------|------------|--------------------------|---------|------------------------|
| 35204/35206 | Black | Soft/Rigid | L/R | X-Small | < 11 cm |
| 35204/35206 | Black | Soft/Rigid | L/R | Small | 11 - 14 cm |
| 35204/35206 | Black | Soft/Rigid | L/R | Medium | 14 - 17 cm |
| 35204/35206 | Black | Soft/Rigid | L/R | Large | 17 - 20 cm |
| 35204/35206 | Svart | Soft | L/R | X-Large | 20 - 23 cm |
| 35205 | Beige | Soft | L/R | X-Small | < 11 cm |
| 35205 | Beige | Soft | L/R | Small | 11 - 14 cm |
| 35205 | Beige | Soft | L/R | Medium | 14 - 17 cm |
| 35205 | Beige | Soft | L/R | Large | 17 - 20 cm |
| 35205 | Beige | Soft | L/R | X-Large | 20 - 23 cm |

Open Wrist

SELECTION® Open Wrist

Selection Open Wrist is a stable yet simple wrist brace in our comfortable Selection-material. The brace is easy to put on with only one hand.

It is made of breathable material and have a shaped plastic spoon on volar side.

The brace has narrow plastic reinforcement on dorsal side and a stabilizing wrap-around-strap.

| Art. nr. | Colour | Left side/Right side | Size | Wrist circumference |
|----------|--------|----------------------|---------|------------------------|
| 35207 | Black | L/R | X-Small | < 11 cm |
| 35207 | Black | L/R | Small | 11 - 14 cm |
| 35207 | Black | L/R | Medium | 14 - 17 cm |
| 35207 | Black | L/R | Large | 17 - 20 cm |
| 35207 | Black | L/R | X-Large | 20 - 23 cm |

SELECTION® Wrist for children

Soft and comfortable wrist support for children. The orthosis is stable with an aluminium-spoon inside and two straps for an easy and stepless closure.

The Thumb straps are adjustable and the edges are soft.

The orthosis has colorful details, chosen to attract the child.

| Art.no | Colour | Left side/ Right side | Size | Wrist circumference | Inside palm | Lenght |
|--------|--------|--------------------------|----------|------------------------|-------------|---------|
| 35305 | Black | L/R | XX-Small | 8,5 - 11 cm | 5-6 cm | 11 cm |
| 35305 | Black | L/R | X-Small | 10 - 12,5 cm | 6-7 cm | 13,5 cm |
| 35305 | Black | L/R | Small | 11,5 - 14 cm | 7-8 cm | 14,5 cm |
| 35305 | Black | L/R | Medium | 12 - 15,5 cm | 8-9 cm | 17 cm |
| 35305 | Black | L/R | Large | 14 - 16,5 cm | 8-9 cm | 17,5 cm |

SELECTION® thumb for children

An soft and comfortable thumb orthosis, intended for children.

The orthosis purpose is to position the hand and wrist in a functional position. It eases pain that comes with reumathical and inflammatorial disorders.

The orthosis can also be used for reducing spasticity within neurological disorders in hand and wrist. The thumb strap is adjustable for optimal fit.

| Art.no | Colour | Left side/ Right side | Size | Wrist Circumference | Inside palm |
|--------|--------|--------------------------|----------|------------------------|-------------|
| 35300 | Black | L/R | XX-Small | 8,5 - 11,5 cm | 5 - 6 cm |
| 35300 | Black | L/R | X-Small | 10 - 12,5 cm | 6 - 7 cm |
| 35300 | Black | L/R | Small | 11,5 - 14 cm | 7 - 8 cm |
| 35300 | Black | L/R | Medium | 12 - 15,5 cm | 8 - 9 cm |
| 35300 | Black | L/R | Large | 14 - 16,5 cm | 8 - 9 cm |



Thumb- ortosis

SELECTION® Thumb Soft & Rigid


Selection Thumb orthosis has been developed from the patterns of our two most popular products. Over the years, these patterns have proven to offer superior fit and function for comfortable thumb stabilization. SELECTION Soft is elastic all the way around and the design keeps the thumb in a relaxed and painfree position. The Rigid version gives more stability due to non-elastic fabric around the thumb section.



| Art. no | Colour | Model | Left side/ Right side | Size | Wrist Circumference |
|---------|--------|-------|--------------------------|----------|---------------------|
| 35200 | Black | Soft | L/R | X-Small | < 13 cm |
| 35200 | Black | Soft | L/R | Small | 13 - 15 cm |
| 35200 | Black | Soft | L/R | Medium | 15 - 17 cm |
| 35200 | Black | Soft | L/R | Large | 17 - 20 cm |
| 35200 | Black | Soft | L/R | X-Large | > 20 cm |
| 35201 | Beige | Soft | L/R | X-Small | <13 cm |
| 35201 | Beige | Soft | L/R | Small | 13 - 15 cm |
| 35201 | Beige | Soft | L/R | Medium | 15 - 17 cm |
| 35201 | Beige | Soft | L/R | Large | 17 - 20 cm |
| 35201 | Beige | Soft | L/R | X-Large | > 20 cm |
| 35202 | Black | Rigid | L/R | X -Small | < 13 cm |
| 35202 | Black | Rigid | L/R | Small | 13 - 15 cm |
| 35202 | Black | Rigid | L/R | Medium | 15 - 17 cm |
| 35202 | Black | Rigid | L/R | Large | 17 - 20 cm |
| 35202 | Black | Rigid | L/R | X-Large | > 20 cm |
| 35203 | Beige | Rigid | L/R | X-Small | < 13 cm |
| 35203 | Beige | Rigid | L/R | Small | 13 - 15 cm |
| 35203 | Beige | Rigid | L/R | Medium | 15 - 17 cm |
| 35203 | Beige | Rigid | L/R | Large | 17 - 20 cm |
| 35203 | Beige | Rigid | L/R | X-Large | > 20 cm |

SELECTION[®]

Accommodative function



When working with the SELECTION line of products, we have chosen a new approach to material selection. New, fresh and with a modern feel it is a combination that meets the users requirements. The SELECTION-material is thin and pliable, let's the skin breath and keeps it cool.

All hook and loop material has been integrated into the material to make the brace as low profile as possible.

This material does not require any edging tape which also contributes to the low profile, but also allows the practitioner to cut away material if neccessary without fraying

