

SELECTION[®]

ULNAR DEVIATION SPLINT

Improve finger alignment to enhance hand function and functional grip.



- Low profile design
- Meets the shape and contours of the hand
- Breathable fabric
- Fabric can be trimmed
- Easy fastening
- Can be worn with Selection Wrist Splint
- Machine washable

SELECTION®

ULNAR DEVIATION



SELECTION® Ulnar Deviation Splint is intended to improve finger alignment and enhance hand function and functional grip.

This hand splint is exclusively designed and produced using our unique three layer Selection fabric. With its low profile design, it meets the contours of the hand, without impeding function. The fabric can be trimmed without fraying, and therefore allows further customization if required. The easy fastening surface allows the hook to be placed anywhere on the splint to optimise the directional pull needed.

The Ulnar Deviation splint supports the arches of the hand whilst providing re-alignment, stabilisation and can contribute to the management of pain. The Orthosis has a slight wider fit ulnarly at the 5th MCP, to counteract the pulling direction in best possible way.

Should you need to provide wrist stabilisation for your patient, this splint can also be worn in combination with a Selection wrist splint.



Intended use

For patients with inflammatory and degenerative diseases, such as rheumatoid arthritis, that has resulted in ulnar deviation. Selection Ulnar Deviation can also be used for hypermobile joints and weak muscles to provide stability. Not intended for use with fixed contractures in the hand.



Item No. Left	Item No. Right	Size	Circumference MCP 2-5
35218 1010	35218 2010	X-Small	14 - 17 cm
35218 1011	35218 2011	Small	17 - 20 cm
35218 1012	35218 2012	Medium	20 - 23 cm
35218 1013	35218 2013	Large	23 - 27 cm