

CROSS™ Hemi Shoulder

CROSS™ Hemi Shoulder is a comfortable and stable textile orthosis, designed to manage subluxation and shoulder pain (see reference literature under Recommended Literature below).

The orthosis is engineered to provide controlled correction of the arm and shoulder position through external rotation and elevation, aimed at alleviating pain, improving posture, and promoting functional movement. In this way, it supports the user's active rehabilitation process.

Recommended Range Of Application

Subluxation, shoulder pain, and dysfunction following conditions such as:
Stroke
Paretic shoulder/arm due to hemiplegia
Brachial plexus injuries

Contraindications

Traumatic shoulder dislocation. Use of Pacemaker. Use of defibrillator (ICD).

Measurements

Measurement 1: Upper arm circumference
Measurement 2: Forearm circumference



Item No.	Model	Size	L/R	Measurement 1	Measurement 2
287171010	CROSS Hemi Shoulder	X-Small	Left	24 - 28 cm	21 - 24 cm
287171011	CROSS Hemi Shoulder	Small	Left	27 - 31 cm	23 - 26 cm
287171012	CROSS Hemi Shoulder	Medium	Left	30 - 35 cm	25 - 29 cm
287171013	CROSS Hemi Shoulder	Large	Left	34 - 39 cm	28 - 32 cm
287171014	CROSS Hemi Shoulder	X-Large	Left	38 - 43 cm	31 - 35 cm
287172010	CROSS Hemi Shoulder	X-Small	Right	24 - 28 cm	21 - 24 cm
287172011	CROSS Hemi Shoulder	Small	Right	27 - 31 cm	23 - 26 cm
287172012	CROSS Hemi Shoulder	Medium	Right	30 - 35 cm	25 - 29 cm
287172013	CROSS Hemi Shoulder	Large	Right	34 - 39 cm	28 - 32 cm
287172014	CROSS Hemi Shoulder	X-Large	Right	38 - 43 cm	31 - 35 cm
287180000	CROSS Hemi Shoulder, Strap	One Size	-	-	-