

Supinator™ PTT Stabiliser

The Supinator is a comfortable and relatively low-profile brace for the treatment of posterior tibial tendon dysfunction (PTTD). With its patent-pending design, the Supinator reduces over-pronation by applying a supinating force. The arch strap works in concert with a repositionable arch pad to lift the medial arch and the heel strap restricts eversion of the calcaneus (heel). The arch pad may be removed if desired to allow the Supinator PTTD brace to be worn with an orthotic. In addition, the Supinator features the Speed Lacer II closure system for ease of application and equal tension across all laces.

Recommended Range Of Application

Posterior tibial tendon dysfunction (PTTD), acquired adult flatfoot deformity
posterior tibial tendon insufficiency.

Contraindications

Lateral ankle sprain, lateral ankle instability.

Measurements

Ankle circumference



Item No.	Model	Size	L/R	Ankle Circumference
281961010	Supinator PTT stabiliser	X-Small	Left	25,4–27,9 cm / 10"-11"
281961011	Supinator PTT stabiliser	Small	Left	27,9–30,5 cm / 11"-12"
281961012	Supinator PTT stabiliser	Medium	Left	30,5–33 cm / 12"-13"
281961013	Supinator PTT stabiliser	Large	Left	33–35,6 cm / 13"-14"
281961014	Supinator PTT stabiliser	X-Large	Left	35,6–38,1 cm / 14"-15"
281962010	Supinator PTT stabiliser	X-Small	Right	25,4–27,9 cm / 10"-11"
281962011	Supinator PTT stabiliser	Small	Right	27,9–30,5 cm / 11"-12"
281962012	Supinator PTT stabiliser	Medium	Right	30,5–33,0 cm / 12"-13"
281962013	Supinator PTT stabiliser	Large	Right	33,0–35,6 cm / 13"-14"
281962014	Supinator PTT stabiliser	X-Large	Right	35,6–38,1 cm / 14"-15"