

# Cubital Tunnel Brace

The Cubital Tunnel Brace is a low-profile night splint that comfortably holds the elbow in a neutral position for the treatment of cubital tunnel syndrome. Ventral design, with no exposed rigid materials, enhances patient compliance during sleep. This cubital tunnel splint features an adjustable aluminum stay that allows the elbow to be positioned at the optimum angle. Unique strap retention system allows for easy application by holding the proximal and distal closure straps within their plastic buckles. The middle strap may be removed if desired for various treatment protocols.



### Recommended Range Of Application

Cubital tunnel syndrome. Ulnar nerve entrapment at the elbow.

### Measurements

Circumference of upper arm at mid bicep

Item No.	Model	L/R	Size	Circumference
285820011	Cubital Tunnel Brace	Bilateral	Small/Medium	15,2–38,1 cm / 6"-15"
285820013	Cubital Tunnel Brace	Bilateral	Large/X-Large	25,4 cm–48,3 cm / 10"-19"