

Elements Body - Vest & Body Suit

Elements Body is an orthosis made of Lycra fabric, composed of various panels. It is designed to give constant proprioceptive feedback through pressure and resistance. The function of Elements Body is to improve posture and stability. This is achieved by providing increased pressure to the body through a snug custom fit. By increasing the deep sensibility, proprioceptive feedback is increased. The proprioceptive feedback may also improve the patient's fine motor skills and functional performance. Various reinforcements can be selected to optimise flexion, extension, pronation, and supination positions.



Recommended Range Of Application

Cerebral Palsy: hypotonia, spasticity, ataxia

MS, Stroke & Acquired Brain Injury: spasticity, ataxia

Muscular Dystrophy, Acquired Spinal Cord Injury and Down's Syndrome: hypotonia

Ataxia: spastic ataxia

Measurements

Measurements are taken ensuring they reflect a snug fit. The order is placed and Elements Body is custom fabricated. The lycra vest can be designed with reinforcements and fastenings as optional extras. The design is determined by the prescribing clinician. Elements Body is available in different colour combinations. Basic materials, reinforcements and stitching can all be chosen at the user's preference.

Item No.	Model
393000000	Elements Body - Vest Without Sleeves
393060012	Elements Body - Body with ¾ sleeve
393010011	Elements Body - Vest With Short Sleeves
393010012	Elements Body - Vest with ¾ sleeve
393010013	Elements Body - Vest With Long Sleeves
393050000	Elements Body - Body Suit Without Sleeves
393060011	Elements Body - Body Suit With Short Sleeves
393060013	Elements Body - Body Suit With Long Sleeves
393500000	Extra Reinforcement Small
393510000	Extra Reinforcement Large
393600000	First Free of Charge Alteration
393610000	Small Alteration
393630000	Large Alteration