

Elements Body – Suit

Elements Body Suit is a textile orthosis that provides continuous proprioceptive feedback across multiple parts of the body. The Elements Body Suit may be particularly suitable when coordination between the trunk, arms, and legs requires support to improve posture, stability, and coordination.

The Elements Body orthosis is made of Lycra elastic fabric and constructed from multiple panels in a design that provides pressure and resistance to enhance posture and stability. By increasing deep sensory input, also referred to as proprioception, fine motor skills as well as functional performance may also be improved.

The suit is ordered and manufactured according to individual measurements and is available with short or long sleeves and legs. Reinforcements and zippers can be added as needed for functional and usability needs. The design and configuration are determined by the prescribing clinician.

Advised minimum wearing time for Elements Body is 6 hours a day. At the beginning of treatment, Elements Body should be worn for shorter periods of time, and gradually increased until the goal of 6 hours per day is reached. The orthosis should not be used for more than 8–10 hours per day. Elements Body should not be worn through the night when patient is sleeping.

Recommended Range Of Application

Cerebral Palsy: hypotonia, spasticity, ataxia

MS, Stroke & Acquired Brain Injury: spasticity, ataxia

Muscular Dystrophy, Acquired Spinal Cord Injury and Down's Syndrome: hypotonia

Ataxia: spastic ataxia



Item No.	Model
393100011	Suit Without Sleeves & Short Legs
393120008	¾-sleeve suit with short legs
393120009	Suit with ¾ sleeves and ¾ length legs
393120010	Suit with ¾ sleeves/long legs
393100013	Suit Without Sleeves & Long Legs
393100012	Suit Without Sleeves / ¾ Legs
393110011	Suit With Sleeves & Short Legs
393110012	Suit With Short Sleeves / ¾ Legs
393110013	Suit With Short Sleeves & Long Legs
393120011	Suit With Long Sleeves & Short Legs
393120013	Suit With Long Sleeves & Long Legs
393120012	Suit With Long Sleeves / ¾ Legs
393500000	Extra Reinforcement Small
393510000	Extra Reinforcement Large
393600000	First Free of Charge Alteration
393610000	Small Alteration
393630000	Large Alteration