

Elements Body - Glove / Sleeve

Elements Body is an orthosis made of Lycra fabric, composed of various panels. It is designed to give constant proprioceptive feedback through pressure and resistance. The function of Elements Body is to improve posture and stability. This is achieved by providing increased pressure to the body through a snug custom fit. By increasing the deep sensibility, proprioceptive feedback is increased. The proprioceptive feedback may also improve the patient's fine motor skills and functional performance. Various reinforcements can be selected to optimise flexion, extension, pronation, and supination positions.

Recommended Range Of Application

Cerebral Palsy: hypotonia, spasticity, ataxia
MS, Stroke & Acquired Brain Injury: spasticity, ataxia
Muscular Dystrophy, Acquired Spinal Cord Injury and Down's Syndrome: hypotonia
Ataxia: spastic ataxia

Measurements

Measurements are taken ensuring they reflect a snug fit. The order is placed and Elements Body is custom fabricated. The lycra vest can be designed with reinforcements and fastenings as optional extras. The design is determined by the prescribing clinician. Elements Body is available in different colour combinations. Basic materials, reinforcements and stitching can all be chosen at the user's preference.



Item No.	Model
393300000	Elements Body - Short Glove Without Fingers
393310000	Elements Body - Short Glove With Fingers
393350000	Elements Body - Long Glove Without Fingers
393360000	Elements Body - Long Glove With Fingers
393400000	Elements Body - Sleeve
393500000	Extra Reinforcement Small
393510000	Extra Reinforcement Large
393600000	First Free of Charge Alteration
393610000	Small Alteration
393630000	Large Alteration