

SELECTION® Childrens Thumb

Selection Childrens Thumb is elastic all the way around and the design keeps the thumb in a relaxed and pain free position. It is designed to guide the thumb and hand into a more functional position to facilitate the gripping ability.

Recommended Range Of Application

Brachial plexus birth palsy (BPBP), Spastic hand in CP children, Arthrogryposis multiplex congenital children, Rheumatic problems.

Contraindications

None known.

Measurements

Measure 1: Wrist circumference

Measure 2: Palm height

SELECTION®



Item No.	Size	L/R	Color	Measurement 1	Measurement 2
353001009	XX-Small	Left	Black	8,5 - 11,5 cm	5 - 6 cm
353001010	X-Small	Left	Black	10 - 12,5 cm	6 - 7 cm
353001011	Small	Left	Black	11,5 - 14 cm	7 - 8 cm
353001012	Medium	Left	Black	12 - 15,5 cm	8 - 9 cm
353001013	Large	Left	Black	14 - 16,5 cm	8 - 9 cm
353002009	XX-Small	Right	Black	8,5 - 11,5 cm	5 - 6 cm
353002010	X-Small	Right	Black	10 - 12,5 cm	6 - 7 cm
353002011	Small	Right	Black	11,5 - 14 cm	7 - 8 cm
353002012	Medium	Right	Black	12 - 15,5 cm	8 - 9 cm
353002013	Large	Right	Black	14 - 16,5 cm	8 - 9 cm