

SELECTION® Wrist Soft

Selection Wrist Soft is a standard wrist brace in a "glove" version with two closing tabs attached to lacing on the ulnar side of the wrist. It has an adjustable thumb strap, and the shape is low cut in the palm area to allow for the bending zone. This model is also well cut around the thumb muscles to allow full movement of the thumb. It comes with a pre-shaped plastic spoon on the palmar side and a soft narrow plastic spoon on the dorsal side. It is intended to be used when the wrist needs support and off-loading in cases such as rheumatic disorders, soft tissue injuries or attritional wear. This model is 18 cm long on the volar side.

Recommended Range Of Application

Inflammations, work related strain or when rheumatic and/or degenerative change are present.

Contraindications

None known.

Measurements

Measure 1 = Wrist circumference

Measure 2 = Mid arm. Circumference 8 cm above the wrist

SELECTION®



Item No.	Model	Size	L/R	Color	Measurement 1	Measurement 2
352041010	Soft	X-Small	Left	Black	< 11 cm	-
352041011	Soft	Small	Left	Black	11 - 14 cm	-
352041012	Soft	Medium	Left	Black	14 - 17 cm	-
352041013	Soft	Large	Left	Black	17 - 20 cm	-
352041014	Soft	X-Large	Left	Black	20 - 23 cm	-
352041015	Soft	XX-Large	Left	Black	23 - 26 cm	-
352041115	Soft	Small Plus	Left	Black	11 - 14 cm	18 - 21 cm
352041125	Soft	Medium Plus	Left	Black	14 - 17 cm	21 - 24 cm
352041135	Soft	Large Plus	Left	Black	17 - 20 cm	24 - 27 cm
352041145	Soft	X-Large Plus	Left	Black	20 - 23 cm	27 - 30 cm
352041155	Soft	XX-Large Plus	Left	Black	23 - 26 cm	30 - 33 cm
352042010	Soft	X-Small	Right	Black	< 11 cm	-
352042011	Soft	Small	Right	Black	11 - 14 cm	-
352042012	Soft	Medium	Right	Black	14 - 17 cm	-
352042013	Soft	Large	Right	Black	17 - 20 cm	-
352042014	Soft	X-Large	Right	Black	20 - 23 cm	-
352042015	Soft	XX-Large	Right	Black	23 - 26 cm	-
352042115	Soft	Small Plus	Right	Black	11 - 14 cm	18 - 21 cm
352042125	Soft	Medium Plus	Right	Black	14 - 17 cm	21 - 24 cm
352042135	Soft	Large Plus	Right	Black	17 - 20 cm	24 - 27 cm
352042145	Soft	X-Large Plus	Right	Black	20 - 23 cm	27 - 30 cm
352042155	Soft	XX-Large Plus	Right	Black	23 - 26 cm	30 - 33 cm
352051010	Soft	X-Small	Left	Beige	< 11 cm	-
352051011	Soft	Small	Left	Beige	11 - 14 cm	-
352051012	Soft	Medium	Left	Beige	14 - 17 cm	-
352051013	Soft	Large	Left	Beige	17 - 20 cm	-
352051014	Soft	X-Large	Left	Beige	20 - 23 cm	-
352051015	Soft	XX-Large	Left	Beige	23 - 26 cm	-
352051115	Soft	Small Plus	Left	Beige	11 - 14 cm	18 - 21 cm
352051125	Soft	Medium Plus	Left	Beige	14 - 17 cm	21 - 24 cm
352051135	Soft	Large Plus	Left	Beige	17 - 20 cm	24 - 27 cm
352051145	Soft	X-Large Plus	Left	Beige	20 - 23 cm	27 - 30 cm
352051155	Soft	XX-Large Plus	Left	Beige	23 - 26 cm	30 - 33 cm
352052010	Soft	X-Small	Right	Beige	< 11 cm	-
352052011	Soft	Small	Right	Beige	11 - 14 cm	-
352052012	Soft	Medium	Right	Beige	14 - 17 cm	-
352052013	Soft	Large	Right	Beige	17 - 20 cm	-
352052014	Soft	X-Large	Right	Beige	20 - 23 cm	-
352052015	Soft	XX-Large	Right	Beige	23 - 26 cm	-
352052115	Soft	Small Plus	Right	Beige	11 - 14 cm	18 - 21 cm
352052125	Soft	Medium Plus	Right	Beige	14 - 17 cm	21 - 24 cm
352052135	Soft	Large Plus	Right	Beige	17 - 20 cm	24 - 27 cm
352052145	Soft	X-Large Plus	Right	Beige	20 - 23 cm	27 - 30 cm
352052155	Soft	XX-Large Plus	-	Beige	23 - 26 cm	30 - 33 cm
352100000	Extra volar spoon	One size	Bilateral	Blue	22 - 168 mm	-