

SELECTION® Wrist Rigid

Selection Wrist Rigid has the same features as the soft version, but also a thumb section added for increased thumb support and stability. The thumb strap can be individually adjusted. It comes with a pre-shaped plastic spoon on the palmar side and a soft narrow plastic spoon on the dorsal side. This product is intended to be used when the first CMC-joint needs support and off-loading in cases such as rheumatic disorders, soft tissue injuries or attritional wear.

Recommended Range Of Application

Inflammations, work related strain or when rheumatic and/or degenerative changes are present.

Contraindications

None known.

Measurements

Wrist circumference

SELECTION®



Item No.	Model	Size	Measurement	L/R	Color
352061010	Rigid	X-Small	< 11 cm	Left	Black
352061011	Rigid	Small	11 - 14 cm	Left	Black
352061012	Rigid	Medium	14 - 17 cm	Left	Black
352061013	Rigid	Large	17 - 20 cm	Left	Black
352061014	Rigid	X-Large	20 - 23 cm	Left	Black
352061015	Rigid	XX-Large	23 - 26 cm	Left	Black
352062010	Rigid	X-Small	< 11 cm	Right	Black
352062011	Rigid	Small	11 - 14 cm	Right	Black
352062012	Rigid	Medium	14 - 17 cm	Right	Black
352062013	Rigid	Large	17 - 20 cm	Right	Black
352062014	Rigid	X-Large	20 - 23 cm	Right	Black
352062015	Rigid	XX-Large	23 - 26 cm	Right	Black