

# SELECTION® Childrens Wrist

SELECTION®



Selection Childrens Wrist is a standard wrist brace in a "glove" version with two closing tabs attached to lacing on the ulnar side of the wrist. It has an adjustable thumb strap, and the shape is low cut in the palm area to allow for the bending zone. This model is also well cut around the thumb muscles to allow full movement of the thumb. It comes with a simple aluminium stay. It is intended to help guide the hand and wrist into a functional position. To reduce pain in combination with inflammatory and rheumatic conditions. To reduce spasticity in hand and wrist in neurological conditions.

### Recommended Range Of Application

Rheumatic conditions. Neurological disorders where wrist position needs mild correction or support. After casting for short periods of use. Inflammatory conditions.

### Contraindications

None known.

### Measurements

Measure 1: Wrist Circumference

Measure 2: Palm height

Item No.	Size	L/R	Color	Length	Measurement 1	Measurement 2
353051009	XX-Small	Left	Black	11 cm	8,5 - 11 cm	5-6 cm
353051010	X-Small	Left	Black	13,5 cm	10 - 12,5 cm	6-7 cm
353051011	Small	Left	Black	14,5 cm	11,5 - 14 cm	7-8 cm
353051012	Medium	Left	Black	17 cm	12 - 15,5 cm	8-9 cm
353051013	Large	Left	Black	17,5 cm	14 - 16,5 cm	8-9 cm
353052009	XX-Small	Right	Black	11 cm	8,5 - 11 cm	5-6 cm
353052010	X-Small	Right	Black	13,5 cm	10 - 12,5 cm	6-7 cm
353052011	Small	Right	Black	14,5 cm	11,5 - 14 cm	7-8 cm
353052012	Medium	Right	Black	17 cm	12 - 15,5 cm	8-9 cm
353052013	Large	Right	Black	17,5 cm	14 - 16,5 cm	8-9 cm