

Ryno Lacer® II wrist orthosis

Highly stable wrist and thumb orthosis with mouldable aluminium splints on the dorsal, volar, and radial sides. The radial thumb splint extends along the full length of the splint and effectively immobilises the thumb. Made from durable suede-like material with a breathable lining. The lacing locking system is easy to open and close with one hand. Ryno Lacer II is an updated and improved version of the Ryno Lacer, featuring enhanced palm padding, increased thumb cushioning, and an updated Velcro system on the underside of the orthosis.



Recommended Range Of Application

Post-traumatic, post-operative, rheumatologic conditions, and spasticity.

Measurements

Wrist circumference

Item No.	Model	Size	L/R	Color	Measurement
285832014	Ryno Lacer II wrist orthosis, Short	X-Large	Right	Black	-
285831014	Ryno Lacer II wrist orthosis, Short	X-Large	Left	Black	-
285832013	Ryno Lacer II wrist orthosis, Short	Large	Right	Black	-
285832012	Ryno Lacer II wrist orthosis, Short	Medium	Right	Black	-
285832011	Ryno Lacer II wrist orthosis, Short	Small	Right	Black	-
285832010	Ryno Lacer II wrist orthosis, Short	X-Small	Right	Black	-
285831013	Ryno Lacer II wrist orthosis, Short	Large	Left	Black	-
285831012	Ryno Lacer II wrist orthosis, Short	Medium	Left	Black	-
285831011	Ryno Lacer II wrist orthosis, Short	Small	Left	Black	-
285831010	Ryno Lacer II wrist orthosis, Short	X-Small	Left	Black	-
285761010	-	X-Small	Left	Black	<14 cm
285761011	-	Small	Left	Black	14-16 cm
285761012	-	Medium	Left	Black	16-18 cm
285761013	-	Large	Left	Black	18-20 cm
285761014	-	X-Large	Left	Black	20-22 cm
285762010	-	X-Small	Right	Black	<14 cm
285762011	-	Small	Right	Black	14-16 cm
285762012	-	Medium	Right	Black	16-18 cm
285762013	-	Large	Right	Black	18-20 cm
285762014	-	X-Large	Right	Black	20-22 cm