

ASO® Ankle Stabiliser

Stable low-profile ankle support with stabilising front lacing. The rigid nylon straps are positioned to mimic taping and prevent the foot from supinating or pronating. The straps can be adjusted depending on the nature of the injury. For added comfort, the orthosis features new reinforced CoolFlex™ padding along the inside, next to the Achilles Tendon.

Recommended Range Of Application

Acute or chronic ankle instability, and for preventive use during sports activities.



Item No.	Model	Size	Measurement	L/R	Color
281860010	ASO Ankle Stabiliser w/Stays	X-Small	25-27,5 cm	Bilateral	Black
281860011	ASO Ankle Stabiliser w/Stays	Small	27,5-30 cm	Bilateral	Black
281860012	ASO Ankle Stabiliser w/Stays	Medium	30-32,5 cm	Bilateral	Black
281860013	ASO Ankle Stabiliser w/Stays	Large	32,5-35 cm	Bilateral	Black
281860014	ASO Ankle Stabiliser w/Stays	X-Large	35-37,5 cm	Bilateral	Black
281860015	ASO Ankle Stabiliser w/Stays	XX-Large	37,5-40 cm	Bilateral	Black
281880009	ASO Ankle Stabiliser w/o Stays	XX-Small	22,5-25 cm	Bilateral	Black
281880010	ASO Ankle Stabiliser w/o Stays	X-Small	25-27,5 cm	Bilateral	Black
281880011	ASO Ankle Stabiliser w/o Stays	Small	27,5-30 cm	Bilateral	Black
281880012	ASO Ankle Stabiliser w/o Stays	Medium	30-32,5 cm	Bilateral	Black
281880013	ASO Ankle Stabiliser w/o Stays	Large	32,5-35 cm	Bilateral	Black
281880014	ASO Ankle Stabiliser w/o Stays	X-Large	35-37,5 cm	Bilateral	Black
281880015	ASO Ankle Stabiliser w/o Stays	XX-Large	37,5-40 cm	Bilateral	Black
281880016	ASO Ankle Stabiliser w/o Stays	XXX-Large	>40 cm	Bilateral	Black