

# NAVIGAIT™

NAVIGAIT™ is an alternative for users with a mild foot drop. The amount of dorsiflexion is easily adjusted with an elastic coil to accommodate the user's specific needs. NAVIGAIT 4-Foot is a sock that makes it possible to use NAVIGAIT without shoes or shoes without an anchor point (e.g., sandals). It is a plastic insert for increased stabilisation and 'stay-up' silicone on the inside to keep it in place.

## Recommended Range Of Application

Intended to support a foot when the ability to actively dorsiflex is reduced due to conditions such as Stroke, Multiple Sclerosis, Post Polio Syndrome.

## Contraindications

NAVIGAT should not be used by patients who present with spasticity or severe oedema.

## Measurements

Length of plastic section



Item No.	Model	Size	Measurement	L/R
284100010	NAVIGAIT	X-Small	23 cm	Bilateral
284100011	NAVIGAIT	Small/Medium	26,5 cm	Bilateral
284100013	NAVIGAIT	Large/X-Large	30 cm	Bilateral
284150011	4-foot	Small	21-23 cm	Bilateral
284150012	4-foot	Medium	23-26 cm	Bilateral
284150013	4-foot	Large	26-29 cm	Bilateral
284120010	Extra strap w/coil	X-Small	-	Bilateral
284120011	Extra strap w/coil	Small/Medium	-	Bilateral
284120013	Extra strap w/coil	Large/X-Large	-	Bilateral
284130010	Extra Strong/Shorter Strap w/coil	X-Small	-	Bilateral
284130011	Extra Strong/Shorter Strap w/coil	Small/Medium	-	Bilateral
284130013	Extra Strong/Shorter Strap w/coil	Large/X-Large	-	Bilateral
284140010	Extra strap set	X-Small	-	Bilateral
284140011	Extra strap set	Small/Medium	-	Bilateral
284140013	Extra strap set	Large/X-Large	-	Bilateral