

NAVIGAIT™

NAVIGAIT™ is an alternative for users with a mild foot drop. The amount of dorsiflexion is easily adjusted with an elastic coil to accommodate the user's specific needs. NAVIGAIT 4-Foot is a sock that makes it possible to use NAVIGAIT without shoes or shoes without an anchor point (e.g., sandals). It is a plastic insert for increased stabilisation and 'stay-up' silicone on the inside to keep it in place.

Recommended Range Of Application

Intended to support a foot when the ability to actively dorsiflex is reduced due to conditions such as Stroke, Multiple Sclerosis, Post Polio Syndrome.

Contraindications

NAVIGAT should not be used by patients who present with spasticity or severe oedema.

Measurements

Length of plastic section



Item No.	Model	Size	Measurement	L/R
284100010	NAVIGAIT	X-Small	23 cm	Bilateral
284100011	NAVIGAIT	Small/Medium	26,5 cm	Bilateral
284100013	NAVIGAIT	Large/X-Large	30 cm	Bilateral
284150011	4-foot	Small	21-23 cm	Bilateral
284150012	4-foot	Medium	23-26 cm	Bilateral
284150013	4-foot	Large	26-29 cm	Bilateral
284120010	Extra strap w/coil	X-Small	-	Bilateral
284120011	Extra strap w/coil	Small/Medium	-	Bilateral
284120013	Extra strap w/coil	Large/X-Large	-	Bilateral
284130010	Extra Strong/Shorter Strap w/coil	X-Small	-	Bilateral
284130011	Extra Strong/Shorter Strap w/coil	Small/Medium	-	Bilateral
284130013	Extra Strong/Shorter Strap w/coil	Large/X-Large	-	Bilateral
284140010	Extra strap set	X-Small	-	Bilateral
284140011	Extra strap set	Small/Medium	-	Bilateral
284140013	Extra strap set	Large/X-Large	-	Bilateral