

# Ypsilon® Flow ½

Ypsilon® Flow ½ is the AFO of choice for active patients who require less stability and more flexibility. Ypsilon® Flow ½ has a lower heel height (approx. 7 mm) which accommodates gentler contours of shoe insoles and offers more clearance for the forefoot in the shoe toe box. Designed to allow optimum range of motion and dynamic toe-off assist, providing the opportunity for functional or potentially functional muscles, tendons, and ligaments to strengthen and prevent or reduce atrophy. The footplate is engineered to allow the orthosis to adapt to and move with the lower leg with less resistance to ground reaction forces, while the proximal ends of the “Y” provide fixation points, eliminating pressure on the tibia crest. Smooth and thin, the interface pad and strap fit close to the leg, ensuring the brace glides smoothly and is barely noticeable under jeans or trousers. Ypsilon® Flow ½ has a heel height that is 7 mm in all sizes.

## Recommended Range Of Application

Ypsilon® Flow is designed to manage footdrop in conditions such as Stroke, Multiple Sclerosis, Post Polio Syndrome, Muscular Dystrophy, Spinal Cord Injuries, Traumatic Brain Injuries, Guillian-Barre Syndrome, Charcot-Marie-Tooth, Myelomeningocele, Neuropathy or Cerebral Palsy.

## Contraindications

Ypsilon® Flow should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits (e.g. quadriceps spasticity, genu valgum or varum, genu recurvatum), severe spasticity.



Item No.	Model	Size	L/R	Footplate Length	Height
288801011	Ypsilon® Flow½	Small	Left	23 cm	35,5 cm
288801012	Ypsilon® Flow½	Medium	Left	24,5 cm	36,5 cm
288801013	Ypsilon® Flow½	Large	Left	27 cm	37,5 cm
288801014	Ypsilon® Flow½	X-Large	Left	28,5 cm	38,5 cm
288802011	Ypsilon® Flow½	Small	Right	23 cm	35,5 cm
288802012	Ypsilon® Flow½	Medium	Right	24,5 cm	36,5 cm
288802013	Ypsilon® Flow½	Large	Right	27 cm	37,5 cm
288802014	Ypsilon® Flow½	X-Large	Right	28,5 cm	38,5 cm
288840012	Tibia pads	S/M/L	Bilateral	-	-
288840014	Tibia pads	X-Large	Bilateral	-	-
288850000	Strap	One size	Bilateral	-	-
998090013	Microfix Ypsilon Flow	-	-	-	-