

# ToeOFF® 2.0

ToeOFF® 2.0 is the updated version of the Original ToeOFF®. This new model has several new features that have been asked for by customers over time. The shape is more organic, with shorter side extensions and a more rounded shape. The surface is smoother and more scratch resistant. The MikroFix is pre-attached on the inside of the brace to accommodate a faster fitting procedure. It is also available with two different strap systems – With D-ring or without D-ring. Please note that the D-ring strapping system is not currently available in Ireland and the UK. ToeOFF® 2.0 has a heel height between 12–15 mm, depending on size.

## Recommended Range Of Application

ToeOFF® 2.0 is designed to manage footdrop in conditions such as Stroke, Multiple Sclerosis, Post Polio Syndrome, Muscular Dystrophy, Spinal Cord Injuries, Traumatic Brain Injuries, Guillian-Barre Syndrome, Charcot-Marie-Tooth, Myelomeningocele, Neuropathy or Cerebral Palsy. ToeOFF® 2.0 is designed to support gait in conditions such as Posterior Tibialis Tendon Dysfunction (PTTD) and toe amputations. ToeOFF® 2.0 can be used for toe amputations.

## Contraindications

ToeOFF® 2.0 should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits (e.g. quadriceps spasticity, genu valgum or varum, genu recurvatum), severe spasticity.

ToeOFF® 2.0



Item No.	Model	Size	L/R	Shoe Size	Footplate Length	Height
289201010	ToeOFF 2.0 Simple Strap	X-Small	Left	32-35	21 cm	36 cm
289201011	ToeOFF 2.0 Simple Strap	Small	Left	35-38	23 cm	38 cm
289201012	ToeOFF 2.0 Simple Strap	Medium	Left	38-42	24,5 cm	40,5 cm
289201013	ToeOFF 2.0 Simple Strap	Large	Left	42-45	27 cm	43 cm
289201014	ToeOFF 2.0 Simple Strap	X-Large	Left	>45	28,5 cm	45 cm
289202010	ToeOFF 2.0 Simple Strap	X-Small	Right	32-35	21 cm	36 cm
289202011	ToeOFF 2.0 Simple Strap	Small	Right	35-38	23 cm	38 cm
289202012	ToeOFF 2.0 Simple Strap	Medium	Right	38-42	24,5 cm	40,5 cm
289202013	ToeOFF 2.0 Simple Strap	Large	Right	42-45	27 cm	43 cm
289202014	ToeOFF 2.0 Simple Strap	X-Large	Right	>45	28,5 cm	45 cm