

BlueROCKER® 2.0

BlueROCKER® is designed to manage foot-drop in conditions such as stroke, Multiple Sclerosis, Post Polio Syndrome, Muscular Dystrophy, Spinal Cord Injuries, Traumatic Brain Injuries, 'Guillian-Barre' Syndrome, Charcot-Marie-Tooth, Myelomeningocele, Neuropathy or Cerebral Palsy. BlueROCKER® 2.0 has a heel height between 12–15 mm, depending on size.

Recommended Range Of Application

BlueROCKER® is designed to support gait in conditions such as Posterior Tibialis Tendon Dysfunction (PTTD) and toe amputations. BlueROCKER® can also be used for partial foot amputations, most proximal level is Chopart.

Contraindications

BlueROCKER® should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits (e.g. quadriceps spasticity, genu valgum or varum, genu recurvatum), severe spasticity.

BLUEROCKER® 2.0



Item No.	Model	Size	L/R	Shoe Size	Footplate Length	Height
289401010	BlueROCKER 2.0	X-Small	Left	32-35	21 cm	36 cm
289401011	BlueROCKER 2.0	Small	Left	35-38	23 cm	38 cm
289401012	BlueROCKER 2.0	Medium	Left	38-42	24,5 cm	40,5 cm
289401013	BlueROCKER 2.0	Large	Left	42-45	27 cm	43 cm
289401014	BlueROCKER 2.0	X-Large	Left	>45	28,5 cm	45 cm
289402010	BlueROCKER 2.0	X-Small	Right	32-35	21 cm	36 cm
289402011	BlueROCKER 2.0	Small	Right	35-38	23 cm	38 cm
289402012	BlueROCKER 2.0	Medium	Right	38-42	24,5 cm	40,5 cm
289402013	BlueROCKER 2.0	Large	Right	42-45	27 cm	43 cm