

# ToeOFF® Short

ToeOFF® Short is specifically developed to accommodate individuals with shorter stature, slightly broader-shaped tibia, and wider foot at the metatarsal area. ToeOFF® Short has a heel height that is 7 mm in all sizes.

## Recommended Range Of Application

ToeOFF® is designed to manage footdrop in conditions such as Stroke, Multiple Sclerosis, Post Polio Syndrome, Muscular Dystrophy, Spinal Cord Injuries, Traumatic Brain Injuries, Guillian-Barre Syndrome, Charcot-Marie-Tooth, Myelomeningocele, Neuropathy or Cerebral Palsy. ToeOFF® is designed to support gait in conditions such as Posterior Tibialis Tendon Dysfunction (PTTD) and toe amputations. BlueROCKER® can also be used for partial foot amputations, most proximal level is Chopart.

## Contraindications

ToeOFF® should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits (e.g. quadriceps spasticity, genu valgum or varum, genu recurvatum), severe spasticity.

ToeOFF® SHORT



Item No.	Model	Size	L/R	Shoe Size	Footplate Length	Height
289901011	ToeOFF Short, Black	Small	Left	32-35	21,5 cm	29,5 cm
289901012	ToeOFF Short, Black	Medium	Left	35-38	23,5 cm	32 cm
289901013	ToeOFF Short, Black	Large	Left	38-42	25,5 cm	34 cm
289902011	ToeOFF Short, Black	Small	Right	32-35	21,5 cm	29,5 cm
289902012	ToeOFF Short, Black	Medium	Right	35-38	23,5 cm	32 cm
289902013	ToeOFF Short, Black	Large	Right	38-42	25,5 cm	34 cm