

Custom ToeOFF®

It is possible to order a Custom ToeOFF® for patients that require a more intimate fit due to: Calf Atrophy/Hypertrophy, Unique Leg Shapes, Alignment Issues, Dynamic Response Changes. This can be done by scanning, casting, or measuring your patient. The product has multiple varied options. These options can be another heel height, stability, or offload for prominences. All options can be found in Product Information under Documents.

ToeOFF®



Item No.	Model
284600000	Custom ToeOFF®
284580000	Custom SoftKIT
284590000	Custom ComfortKIT