

# Custom BlueROCKER®

It is possible to order a Custom BlueROCKER® for patients that require a more intimate fit due to: Calf Atrophy/Hypertrophy, Unique Leg Shapes, Alignment Issues, Dynamic Response Changes. This can be done by scanning, casting, or measuring your patient. The product has multiple varied options. These options can be another heel height, stability, or offload for prominences. All options are presented in PDF below.

BLUEROCKER®



Item No.	Model
284800000	Custom BlueROCKER®
284580000	Custom SoftKIT
284590000	Custom ComfortKIT