

# ToeOFF® 2.0 Trial Brace

This product is permanently labelled "Trial Brace". It is only to be used for patient assessment, to determine ToeOFF® product and/or size to achieve optimum gait for the patient. We recommend the first assessment to be done in parallel bars. Includes SoftKIT – Applied!

## Recommended Range Of Application

This product is intended only for temporary use.

ToeOFF®



| Item No.  | Model                                | Size    | L/R   | Shoe Size | Footplate Length | Height  |
|-----------|--------------------------------------|---------|-------|-----------|------------------|---------|
| 289211010 | ToeOFF 2.0 Trial Brace Simple Straps | X-Small | Left  | 32-35     | 21 cm            | 36 cm   |
| 289211011 | ToeOFF 2.0 Trial Brace Simple Straps | Small   | Left  | 35-38     | 23 cm            | 38 cm   |
| 289211012 | ToeOFF 2.0 Trial Brace Simple Straps | Medium  | Left  | 38-42     | 24,5 cm          | 40,5 cm |
| 289211013 | ToeOFF 2.0 Trial Brace Simple Straps | Large   | Left  | 42-45     | 27 cm            | 43 cm   |
| 289211014 | ToeOFF 2.0 Trial Brace Simple Straps | X-Large | Left  | >45       | 28,5 cm          | 45 cm   |
| 289212010 | ToeOFF 2.0 Trial Brace Simple Straps | X-Small | Right | 32-35     | 21 cm            | 36 cm   |
| 289212011 | ToeOFF 2.0 Trial Brace Simple Straps | Small   | Right | 35-38     | 23 cm            | 38 cm   |
| 289212012 | ToeOFF 2.0 Trial Brace Simple Straps | Medium  | Right | 38-42     | 24,5 cm          | 40,5 cm |
| 289212013 | ToeOFF 2.0 Trial Brace Simple Straps | Large   | Right | 42-45     | 27 cm            | 43 cm   |
| 289212014 | ToeOFF 2.0 Trial Brace Simple Straps | X-Large | Right | >45       | 28,5 cm          | 45 cm   |