

CoverKIT™ 2.0

CoverKIT™ 2.0 is another new smart interface innovation from Allard. This new interface is designed to fit the 2.0 and 2 ½ models of ToeOFF® and BlueROCKER®. CoverKIT™ 2.0 combines our SoftKIT™ interface with a thin textile front to form a sleeve that slides over the anterior shell. Smooth, sleek appearance. The kit includes two narrow neoprene pads if there is a need to create a pressure-free channel for tibia prominences. If needed, the original straps, supplied with the brace, can be worn together with the CoverKIT™.



Item No.	Model	Size	Color	Calf circumference
287650010	coverKIT™ 2.0	X-Small	Black	31 - 34 cm
287650011	coverKIT™ 2.0	Small	Black	34 - 37 cm
287650012	coverKIT™ 2.0	Medium	Black	37 - 41 cm
287650013	coverKIT™ 2.0	Large	Black	41 - 46 cm
287650014	coverKIT™ 2.0	X-Large	Black	46 - 51 cm