

## SoftKIT™

Every ToeOFF® and BlueROCKER® should have a soft interface between the orthosis and the tibial crest. SoftKIT is a padding system that includes two narrow neoprene pads to create a pressure-free channel to help relieve any potential discomfort from bony prominences. The kit includes pre-cut pressure-sensitive MikroFIX tabs to easily affix the pad to ToeOFF® or BlueROCKER®. An extra pad is included, so one can be worn while the other is laundered.



Item No.	Model	Size	Color
283900010	Interface for ToeOFF & BlueROCKER	X-Small	Black
283900011	Interface for ToeOFF & BlueROCKER	Small	Black
283900012	Interface for ToeOFF & BlueROCKER	Medium	Black
283900013	Interface for ToeOFF & BlueROCKER	Large/X-Large	Black
283960011	MicroFIX Kit for ToeOFF/BlueROCKER	X-Small/Small	Black
283960013	MicroFIX Kit for ToeOFF/BlueROCKER	Medium/Large/X-Large	Black
283950011	Vertical Tibia Pad	X-Small/Small	Black
283950013	Vertical Tibia Pad	Medium/Large/X-Large	Black