

KiddieROCKER®

KiddieROCKER® is developed to meet children with more severe problems, and not only foot drop. The kids may have bilateral weakness and need more proximal stability through Stance Phase. Just as the other Paediatric Allard models, KiddieROCKER® is dynamic graded in stability in the different sizes. This model also gives an increased stability in M/L and A/P plane.

Recommended Range Of Application

KiddieROCKER® is designed to support footdrop, gait deviations secondary to proprioceptive deficit (either unstable or low-tone gait), toe-walking with no midfoot collapse, low tone crouch gait in conditions such as Spina Bifida, Cerebral Palsy, Muscular Dystrophy.

Contraindications

KiddieROCKER® should not be used when patients present with: Lacking ROM towards dorsiflexion (need at least 5o dorsiflexion past neutral). Very rigid foot structure, Quadriceps spasticity, excessive postural Genu Valgum or Genu Varum, fixed/ non correctable postural Pes Valgus or Pes Varus.

kiddieROCKER®



Item No.	Model	Size	L/R	Footplate Length	Height
289021011	KiddieROCKER	Small	Left	16 cm	22 cm
289021012	KiddieROCKER	Medium	Left	18 cm	25,7 cm
289021013	KiddieROCKER	Large	Left	20 cm	29,5 cm
289021014	KiddieROCKER	X-Large	Left	21 cm	31,5 cm
289022011	KiddieROCKER	Small	Right	16 cm	22 cm
289022012	KiddieROCKER	Medium	Right	18 cm	25,7 cm
289022013	KiddieROCKER	Large	Right	20 cm	29,5 cm
289022014	KiddieROCKER	X-Large	Right	21 cm	31,5 cm