

SELECTION® Patella

This brace is especially designed for treatment of different kinds of patella related problems. The patella pad can be moved into different positions – medially, laterally or underneath the patella. The brace gives a light stabilization and has stainless steel spiral stays on the sides.

Recommended Range Of Application

Patella dislocation, patellofemoral pain syndrome, Schlatter, Jumpers Knee, anterior knee pain.

Measurements

- Measure 1: 15 cm below patella
- Measure 2: Mid-patella
- Measure 3: 15 cm above patella

SELECTION®



Item No.	Size	Color	Measurement 1	Measurement 2	Measurement 3
357040010	X-Small	Black	28-31 cm	30-33 cm	34-39 cm
357040011	Small	Black	31-34 cm	33-36 cm	39-43 cm
357040012	Medium	Black	34-37 cm	36-39 cm	43-47 cm
357040013	Large	Black	37-40 cm	39-42 cm	47-51 cm
357040014	X-Large	Black	40-43 cm	42-46 cm	51-56 cm