

# SELECTION® Knee Minor

Knee Minor is a tube version with stainless steel spiral stays on each side. It gives mild support and good compression. Minor is adjustable on the front of the thigh and on the calf. It also has an anatomically shaped patella opening and silicone "Stay-Up" on the inside to prevent migration.

### Recommended Range Of Application

Knee pain in weight bearing, minor injuries such as stree in ligaments or mild sprains.

### Contraindications

Severe knee injuries such as ACL / PCL ruptures. Be vigilant for severe problems with swelling and / or poor circulation.

### Measurements

- Measure 1: 15 cm below patella
- Measure 2: Mid-patella
- Measure 3: 15 cm above patella



Item No.	Size	Color	Measurement 1	Measurement 2	Measurement 3
357000010	X-Small	Black	28-31 cm	30-33 cm	34-39 cm
357000011	Small	Black	31-34 cm	33-36 cm	39-43 cm
357000012	Medium	Black	34-37 cm	36-39 cm	43-47 cm
357000013	Large	Black	37-40 cm	39-42 cm	47-51 cm
357000014	X-Large	Black	40-43 cm	42-46 cm	51-56 cm
357000015	XX-Large	Black	43-47 cm	46-50 cm	56-60 cm