

# SELECTION® Knee Minor

Knee Minor is a tube version with stainless steel spiral stays on each side. It gives mild support and good compression. Minor is adjustable on the front of the thigh and on the calf. It also has an anatomically shaped patella opening and silicone "Stay-Up" on the inside to prevent migration.

## Recommended Range Of Application

Knee pain in weight bearing, minor injuries such as stree in ligaments or mild sprains.

## Contraindications

Severe knee injuries such as ACL / PCL ruptures. Be vigilant for severe problems with swelling and / or poor circulation.

## Measurements

Measure 1: 15 cm below patella

Measure 2: Mid-patella

Measure 3: 15 cm above patella



| Item No.  | Size     | Color | Measurement 1 | Measurement 2 | Measurement 3 |
|-----------|----------|-------|---------------|---------------|---------------|
| 357000010 | X-Small  | Black | 28-31 cm      | 30-33 cm      | 34-39 cm      |
| 357000011 | Small    | Black | 31-34 cm      | 33-36 cm      | 39-43 cm      |
| 357000012 | Medium   | Black | 34-37 cm      | 36-39 cm      | 43-47 cm      |
| 357000013 | Large    | Black | 37-40 cm      | 39-42 cm      | 47-51 cm      |
| 357000014 | X-Large  | Black | 40-43 cm      | 42-46 cm      | 51-56 cm      |
| 357000015 | XX-Large | Black | 43-47 cm      | 46-50 cm      | 56-60 cm      |