

SELECTION® Knee Minor Wrap

SELECTION®Knee Minor Wrap is a "wrap-around-model" with stainless steel spiral stays on each side. It gives mild support and good compression. It has an anatomically shaped patella opening and silicone "Stay-Up" on the inside to prevent migration. Minor Wrap also comes in PLUS sizes to accommodate larger thigh circumferences.

Recommended Range Of Application

Knee pain in weight bearing, minor injuries such as stretch in ligaments or mild sprains.

Contraindications

Severe knee injuries such as ACL / PCL ruptures. Be vigilant for severe problems with swelling and / or poor circulation.

Measurements

- Measure 1: 15 cm below patella
- Measure 2: Mid-patella
- Measure 3: 15 cm above patella



Item No.	Model	Size	Color	Measurement 1	Measurement 2	Measurement 3
357010010	SELECTION® Knee Minor Wrap	X-Small	Black	25-28 cm	27-30 cm	31-36 cm
357010011	SELECTION® Knee Minor Wrap	Small	Black	28-31 cm	30-33 cm	36-40 cm
357010115	SELECTION® Knee Minor Wrap	Small Plus	Black	28-31 cm	30-33 cm	40-44 cm
357010012	SELECTION® Knee Minor Wrap	Medium	Black	31-34 cm	33-36 cm	40-44 cm
357010125	SELECTION® Knee Minor Wrap	Medium Plus	Black	31-34 cm	33-36 cm	44-48 cm
357010013	SELECTION® Knee Minor Wrap	Large	Black	34-37 cm	36-39 cm	44-48 cm
357010135	SELECTION® Knee Minor Wrap	Large Plus	Black	34-37 cm	36-39 cm	48-53 cm
357010014	SELECTION® Knee Minor Wrap	X-Large	Black	37-40 cm	39-43 cm	48-53 cm
357010145	SELECTION® Knee Minor Wrap	X-Large Plus	Black	37-40 cm	39-43 cm	53-60 cm
357010015	SELECTION® Knee Minor Wrap	XX-Large	Black	40-46 cm	43-50 cm	53-60 cm
357010155	SELECTION® Knee Minor Wrap	XX-Large Plus	Black	40-46 cm	43-50 cm	60-67 cm