

# SELECTION® Knee Optima

This model is a "pull-on" version. It combines a range of motion control with medial/lateral stability and compressive support for knee instability. It also gives assistance in reducing hyperextension, adding M-L stability and can assist in A-P instability. Furthermore, it has an anatomically shaped patella opening and silicone "Stay-Up" on the inside to prevent migration.

### Recommended Range Of Application

Mild to Moderate knee hyperextension, Weak ACL, Partial ACL rupture, Mild to moderate M-L instability, Mild Osteoarthritis, Post surgery.

### Contraindications

Severe knee injuries such as ACL/PCL ruptures. Be vigilant for severe problems with swelling and / or poor circulation.

### Measurements

- Measure 1: 15cm below patella
- Measure 2: Mid-patella
- Measure 3: 15cm above patella

SELECTION®



Item No.	Size	Color	Measurement 1	Measurement 2	Measurement 3
357020010	X-Small	Black	28-31 cm	30-33 cm	34-39 cm
357020011	Small	Black	31-34 cm	33-36 cm	39-43 cm
357020012	Medium	Black	34-37 cm	36-39 cm	43-47 cm
357020013	Large	Black	37-40 cm	39-42 cm	47-51 cm
357020014	X-Large	Black	40-43 cm	42-46 cm	51-56 cm
357020015	XX-Large	Black	43-47 cm	46-50 cm	56-60 cm