

SELECTION® Knee Optima Wrap

This model is a "wrap-around-model". It combines a range of motion control with medial/lateral stability and compressive support for knee instability. It also gives assistance in reducing hyperextension, adding M-L stability and can assist in A-P instability. Likewise, it has an anatomically shaped patella opening and silicone "Stay-Up" on the inside to prevent migration. Optima Wrap also comes in PLUS sizes to accommodate larger thigh circumferences.

Recommended Range Of Application

Mild to Moderate knee hyperextension, Weak ACL, Partial ACL rupture, Mild to moderate M-L instability, Mild Osteoarthritis, Post surgery.

Contraindications

Severe knee injuries such as ACL/PCL ruptures. Be vigilant for severe problems with swelling and / or poor circulation.

Measurements

- Measure 1: 15 cm below patella
- Measure 2: Mid-ptaella
- Measure 3: 15 cm above patella

SELECTION®



Item No.	Size	Color	Measurement 1	Measurement 2	Measurement 3
357030010	X-Small	Black	25-28 cm	27-30 cm	31-36 cm
357030011	Small	Black	28-31 cm	30-33 cm	36-40 cm
357030115	Small Plus	Black	28-31 cm	33-33 cm	40-44 cm
357030012	Medium	Black	31-34 cm	33-36 cm	40-44 cm
357030125	Medium Plus	Black	31-34 cm	33-36 cm	44-48 cm
357030013	Large	Black	34-37 cm	36-39 cm	44-48 cm
357030135	Large Plus	Black	34-37 cm	36-39 cm	48-53 cm
357030014	X-Large	Black	37-40 cm	39-43 cm	48-53 cm
357030145	X-Large Plus	Black	37-40 cm	39-43 cm	53-60 cm
357030015	XX-Large	Black	40-46 cm	43-50 cm	53-60 cm
357030155	XX-Large Plus	Black	40-46 cm	43-50 cm	60-67 cm