

FreeMotion Joints

This free, open joint is used in conjunction with MultiMotion Modular and MultiStatic treatment joints to ensure torsional rigidity of the orthosis and provide optimal function. Made of stainless steel.



Item No.	Model	Length	Diameter	Width	Height	Weight
696970000	FreeMotion Free Movable Joint Small	93 mm	18 mm	12 mm	11 mm	40 g
696980000	FreeMotion Free Movable Joint T-Bar Small	97 mm	18 mm	12 mm	11 mm	72 g
696950000	FreeMotion Free Movable Joint Small	400 mm	18 mm	12 mm	11 mm	102 g
696990000	FreeMotion Free Movable Joint T-Bar Small	250 mm	18 mm	12 mm	11 mm	106 g
697000000	FreeMotion Free Movable Joint Medium	97 mm	22 mm	16 mm	11 mm	52 g
697030000	FreeMotion Free Movable Joint T-Bar Medium	105 mm	22 mm	16 mm	11 mm	88 g
696960000	FreeMotion Free Movable Joint Medium	400 mm	24 mm	16 mm	8 mm	114 g
697040000	FreeMotion Free Movable Joint T-Bar Medium	260 mm	22 mm	16 mm	11 mm	134 g